CLEANING AND STERILISATION MANUAL

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

Plus / Flex Line

Copyright

© Mectron S.p.A. 2021. All rights reserved. No part of this document can be reproduced in any form without the written consent of the copyright owner.

EN

SUMMARY

1	Introduction	1
2	Disassembly of Parts	2
3	Cleaning and Sterilisation	5
3.1	Cleaning	6
3.1.1	Cleaning the Device Console and the Foot Pedal	6
3.1.1.1	Preparation	6
3.1.1.2	Materials Required	6
3.1.1.3	Cleaning and Disinfection	7
3.1.2	Pre-Cleaning	8
3.1.2.1	Preparation	8
3.1.2.2	Materials Required	8
3.1.2.3	Immersion in Enzymatic Detergent	9
3.1.3	Manual Cleaning of Accessories	12
3.1.3.1	Preparation	12
3.1.3.2	Materials Required	12
3.1.3.3	Immersion in Enzymatic Detergent	13
3.1.3.4	Immersion in Neutral Detergent	15
3.1.3.5	Ultrasonic Bath	17
3.1.4	Automatic Cleaning of Accessories	18
3.1.4.1	Washer - Disinfector	18
3.2	Cleaning Check	19
3.3	Drying and Lubrication of Accessories	20
3.4	Sterilisation	21
3.4.1	Preparation	21
3.4.2	Materials Required	21

PAGE INTENTIONALLY LEFT BLANK

1 INTRODUCTION

Read this manual carefully before proceeding to install, use, maintain or carry out any other interventions on the device. Make sure it is available at all times in the reprocessing rooms of the device and its accessories.

IMPORTANT: to avoid any injury to people or damage to property, please pay careful attention to all the “safety precautions” included in this manual. Depending on the level of risk involved, the safety requirements are classified as follows:

 **WARNING:** (always related to personal injury)

 **CAUTION:** (relating to possible damage to objects)

A **NOTE** identifies special information that may allow for simpler system maintenance or clarify or emphasise important instructions.

The main aim of this manual is to protect patients, surgeons and hospital personnel from possible infections caused by contaminated instruments. It is therefore essential to pay careful attention when performing each step of the reprocessing procedure. These instructions provide detailed information to allow hospital personnel to safely manage all operations and carry out effective cleaning of the:

- Piezosurgery® Plus/Flex Line device;
- Foot pedal;

and effective cleaning and sterilisation of:

- Handpiece;
- Torque wrench.

The hospital or site that physically carries out the procedures is responsible for cleaning and sterilising the instruments.

The procedures and parameters described in this manual have been validated by laboratory tests carried out by an independent entity. These tests have confirmed that the sterility level (SAL) of 10^{-6} is only achieved if the instructions described in this manual are followed perfectly.

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

 **CAUTION:** The instructions reported in the following are intended for qualified personnel able to treat medical devices for their re-use. These instructions have been validated by Mectron S.p.a. The site remains responsible for ensuring that reprocessing is carried out in such a way as to achieve the desired results, by using the instructions, the equipment and the materials described here, and qualified personnel. The process must be regularly monitored and validated. Any deviation applied by the reprocessing site with respect to that specified in these instructions must be validated and the effectiveness and potential adverse consequences assessed.

NOTE: The following PIEZOSURGERY® Plus/Flex Line accessories are supplied non-sterile and must be cleaned and sterilised at first use and at each re-use: handpieces, torque wrenches.

 **WARNING:** The procedures described in this manual do not apply to sterile, single-use accessories. These accessories cannot be reused. Single-use parts are supplied in sterile packaging and are not intended to be re-used. The “cleaning and sterilisation manual” must be used together with the “use and maintenance manual” supplied with the Piezosurgery® Plus/Flex Line unit.

LIMITATIONS TO REPEATED PROCESSES:

Repetition of processes has a minimal effect on these instruments. The end of the life is generally a result of the wear and tear and damages suffered through use. Wear and tear of instruments must be determined through visual examination.

⚠ WARNING: Cleaning and sterilisation of new or repaired tools. All new and repaired device accessories are supplied in non-sterile conditions. At FIRST USE, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in this manual.

⚠ WARNING: Infection control. Special care must be taken in handling contaminated inserts as they are a cause of a potential danger of infection. Hospital personnel must wear suitable personal protection equipment.

⚠ WARNING: Handle sharp and pointed instruments carefully.

2 DISASSEMBLY OF PARTS

Before proceeding with the cleaning procedures described in *Chapter 3 on page 5*, disconnect all the device accessories and components.

⚠ WARNING: Disconnect the device from the power supply using the GENERAL SWITCH on the device console and disconnect it from the electricity supply. Disconnect all leads from the device (handpiece cord and foot pedal).


⚠ WARNING: Supplementary equipotential cable. If present, disconnect the supplementary equipotential cable before cleaning or sterilising.

Disconnect the device foot pedal: press the foot pedal connector release sleeve and pull the connector backwards;


⚠ CAUTION: Do not attempt to unscrew or rotate the connector when disconnecting: the connector may become damaged.

⚠ CAUTION: When disconnecting the foot pedal cable, always hold the cord connector only. Never pull the cord.

1

2


Disconnect the handpiece from the device;

⚠ **CAUTION:** Grip the connector only when carrying out this operation. Do not carry out this manoeuvre using the cord. Do not attempt to unscrew the connector when disconnecting the handpiece. The connector may become damaged.

⚠ **CAUTION:** The handpiece and its cord cannot be separated.

3


Remove the irrigation kit as follows:

- Remove the spike from the irrigation bag;


NOTE: Close the irrigation kit clamps before disconnecting it from the irrigation bag.

4


- Open the peristaltic pump cover completely and remove the tube;
- Close the peristaltic pump cover;

5


- Remove the tube and the clips from the handpiece cord;

6

- Disconnect the tube from the nozzle on the back of the handpiece;
- Dispose of the irrigation kit in accordance with the current regulations governing hospital waste;

7


8

Remove the insert from the handpiece using the torque wrench supplied with the device;

Dispose of the insert in accordance with current regulations governing hospital waste;

⚠ CAUTION: Only use the Piezosurgery® Plus/Flex Line torque wrench to tighten and remove the insert on/from the handpiece. Do not use any other tools, such as pliers, scissors, etc.


9

Unscrew the insert cover or metallic front cone of the handpiece by twisting in an anticlockwise direction;

Dispose of the insert cover in accordance with current national regulations governing hospital waste;

⚠ CAUTION: The metallic cone should not be disposed of; rather it should be reprocessed with the other Piezosurgery® Plus/Flex Line instruments.


10

FOR INSERTS WITH EXTENSION

Remove the extension from the handpiece using the torque wrench supplied;

Dispose of the extension in accordance with current national regulations governing hospital waste;


Connect the protection cap to the end part of the cord, at the connector cable gland;
Insert the cap into the handpiece cord connector.

⚠ **CAUTION:** Always install the protection cap before cleaning.

⚠ **CAUTION:** Before starting cleaning, always visually inspect the protection cap for any damage. Replace if any anomalies are seen.

11

EN

3 CLEANING AND STERILISATION

This table is indicative, only. For the complete cleaning and sterilization procedures of the individual parts, refer to the chapters indicated in the table.

⚠ **CAUTION:** Methods not reported in the table below must not be used.

⚠ **WARNING: Sterile single-use.** Sterile single-use objects must only be used in one surgical procedure and on one patient. Single-use objects must not be reused. Separate and dispose of each single-use object in accordance with current regulations governing hospital waste.

Device Console and Foot Pedal		
Phase	Chapter	Procedure
I	3.1.1	Manual cleaning with detergent

Accessories (Handpiece and Torque Wrench)				
Phase	Chapter	Procedure	Handpiece	Torque Wrench
II	3.1.2	Pre-Cleaning ^{a)}	Immersion into enzymatic detergent	
III	3.1.3	Manual cleaning	Immersion into enzymatic detergent	
			Immersion into detergent	Ultrasonic washing
	3.1.4	Automatic cleaning	Washer - disinfectant with detergent	
IV	3.2	Cleaning check	Visual inspection	
V	3.3	Drying and lubrication	Drying	Drying and Lubrication
VI	3.4	Sterilisation	Packaging, Sterilization and Storage	

a) Optional.

3.1 Cleaning

Device console and foot pedal:

The only cleaning method applicable is the manual method.

Accessories

(Handpiece and torque wrench):

The operator needs to choose from the following cleaning procedures:

- Manual cleaning;
- Automatic cleaning.

The procedures applicable to the Piezosurgery® Plus/Flex Line accessories described in the following manual have been validated by an independent body.

3.1.1 Cleaning the Device Console and the Foot Pedal

3.1.1.1 Preparation

Check that all the following accessories have been removed/disconnected from the device (see *Chapter 2 on page 2*):

- electric power supply lead;
- irrigation kit;
- handpiece;
- irrigation bag support rod;
- supplementary equipotential cable connector (if present);
- foot pedal.

⚠ WARNING: Always turn the device off by means of the switch and disconnect it from the electric grid before cleaning or disinfecting.

⚠ WARNING: The device console is not protected against the penetration of liquids. Do not spray liquids directly onto the surface of the device.

⚠ CAUTION: Do not immerse the device console or foot pedal in liquids and/or solutions of other types.

⚠ CAUTION: Do not sterilise the device console or the foot pedal.

⚠ CAUTION: Do not use machinery/devices to clean the device console or foot pedal.

⚠ CAUTION: Do not rinse the electric contacts of the foot pedal under running water.

3.1.1.2 Materials Required

- Clean, soft lint-free cloths;
- Detergent (pH 6-9) and, if necessary, a non-aggressive neutral pH (pH 7) disinfectant solution.

⚠ CAUTION: For disinfection, if required, it is recommended to use water-based disinfectant solutions, with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated as they can discolour and/or damage plastic materials. This also applies to chemicals such as acetone and alcohol.

3.1.1.3 Cleaning and Disinfection


Clean the device and foot pedal surfaces with a clean, soft lint-free cloth, dampened with a detergent solution (pH 6-9), following the manufacturer's instructions;

For disinfection, if required, it is recommended to use a non-aggressive neutral pH (pH7) solution prepared according to the manufacturer's instructions;


Dry the device body using a clean, non-abrasive, soft lint-free cloth;

Dry the foot pedal using a clean, non-abrasive, soft lint-free cloth.


1

2

3

3.1.2 Pre-Cleaning

The "pre-cleaning" cycle is optional and preparatory to the chosen cleaning method (manual or automatic) and is recommended in the following cases:

- if the accessories to be cleaned are particularly contaminated and/or dirty;
- if a considerable time has passed since last use and pre-cleaning.

The person in charge of reprocessing activities will assess the need for a possible pre-cleaning cycle according to the conditions of the accessories, tools and specific requirements.

3.1.2.1 Preparation

The process applies to the following accessories:

- Handpiece;
- Metallic front cone;
- Torque wrench.

After having disconnected the accessories, following the instructions provided in *Chapter 2 on page 2* proceed with the cleaning procedure described in the following.

⚠ **CAUTION:** Before cleaning the handpiece, remove the insert.

⚠ **CAUTION:** Always install the protection cap before cleaning.

⚠ **CAUTION:** Before proceeding with the cleaning, always visually inspect the protection cap for any damage. Replace it if any anomalies are detected.

⚠ **CAUTION:** Do not immerse the handpiece in ultrasonic tank.

3.1.2.2 Materials Required

- Enzymatic detergent;
- Tap water;
- Suitably-sized container;
- Clean soft bristled, nylon brush;
- Demineralized/distilled water;
- Syringe.

⚠ **CAUTION:** Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated as they can discolour and/or damage plastic materials. This also applies to chemicals such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.

⚠ **CAUTION: Cleaning instruments.** Use appropriate cleaning instruments as specified in this manual. Never use brushes, metallic sponges, abrasive instruments or sharp objects to remove residues during manual cleaning. These may ruin the product surface and permanently and definitively damage the instruments. Only use brushes with soft, nylon bristles.

⚠ **CAUTION:** When cleaning, do not use organic solvents such as acetone or isopropyl alcohol to clean the instruments. Saline cleaning agents containing aldehydes, mercury, active chlorine, chloride, bromine or bromide, iodine or iodide or fixing or corrosive proteins must never be used. Do not place or soak instruments in a rinsing solution. Do not use solvents, lubricants or other products with chemical characteristics that differ from those specified. Use of these products can cause the instruments to malfunction or the release of undesired substances at the operating site.

⚠ **CAUTION:** Dispose of the enzymatic detergent correctly once used. Do not recycle.

3.1.2.3 Immersion in Enzymatic Detergent

This procedure applies to the following accessories:

- Handpiece;
- Torque wrench.
- Metallic front cone;

Prepare an enzymatic detergent solution

b) according to the manufacturer's instructions;


- b) Process validated by independent body with enzymatic detergent Enzymec (Mectron).

Place the handpiece, metallic front cone and torque wrench in a container. Add sufficient enzymatic detergent solution to completely cover the accessories;


! CAUTION: Do not immerse the handpiece if warm.

Leave the handpiece, metallic front cone and torque wrench soaking for 10 minutes at 40 °C ±2; This will reduce the organic residues on the instruments;


1

2

3


HANDPIECE:


Use a single-use syringe (min. 20 ml.) to inject the enzymatic detergent solution into the cavities and internal channels of the handpiece three times in order to effectively remove the residues from the internal channel surfaces;

While immersed into the enzymatic detergent solution, gently brush the surface of the metallic front cone and wrench, using a clean soft bristle, nylon brush to eliminate all visible traces of soil;

4

5

Using a clean soft bristled, nylon brush under running water (minimum quality: tap water, temperature $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$) clean thoroughly but without damaging the surfaces of the handpiece, front cone and torque wrench. Repeat for at least 30 seconds, until all visible traces of soil have been removed.

Also rinse the handpiece cord;

⚠ WARNING: The following parts must be brushed carefully:

Handpiece

- Threaded part onto which the inserts are screwed;
- Threaded part onto which the metallic front cone is screwed;
- Parts adjacent to threaded parts;

Metallic front cone

- The internal threaded part;

Torque wrench

- External metallic ring;
- Internal channels and cavities.




NOTE: Rinse the cavities, grooves and fissures and critical parts thoroughly.

⚠ CAUTION: To prevent the accumulation of enzymatic detergent residues, rinse the instrument surfaces thoroughly.


Rinse the internal channels, injecting cold distilled/demineralised water from a disposable 20 ml syringe three times, to eliminate all enzymatic detergent residues.

NOTE: Rinse the channel within the handpiece by injecting the liquid from the back, taking care to use a waste water container to avoid sprays. Leave the water to drain through the internal channel. Repeat the procedure until the rinsing water is free from residues.

6

7

3.1.3 Manual Cleaning of Accessories

3.1.3.1 Preparation

The process applies to the following accessories:

- Handpiece;
- Metallic front cone;
- Torque wrench.

After having disconnected the accessories, following the instructions provided in *Chapter 2 on page 2* proceed with the cleaning procedure described in the following.

⚠ **CAUTION:** Before cleaning the handpiece, remove the insert.

⚠ **CAUTION:** Always install the protection cap before cleaning.

⚠ **CAUTION:** Before proceeding with the cleaning, always visually inspect the protection cap for any damage. Replace it if any anomalies are detected.

⚠ **CAUTION:** Do not immerse the handpiece in ultrasonic tank.

3.1.3.2 Materials Required

- Enzymatic detergent;
- Neutral pH cleaner;
- Tap water;
- Suitably-sized container;
- Clean soft bristled, nylon brush;
- Ultrasound tank;
- Demineralised/distilled water;
- Syringe;
- Water jet gun.

⚠ **WARNING:** Run the reprocessing procedures at the end of the intervention. Do not allow contaminated instruments to dry before reprocessing.

⚠ **CAUTION:** During cleaning, do not use detergents with a pH outside the range of 6-9. Do not use organic solvents such as acetone or isopropyl alcohol to clean the instruments. Saline cleaning agents containing aldehydes, mercury, active chlorine, chloride, bromine or bromide, iodine or iodide or fixing or corrosive proteins must never be used. Do not place or soak instruments in a rinsing solution. Do not use solvents, lubricants or other products with chemical characteristics that differ from those specified. Use of these products can cause the instruments to malfunction or the release of undesired substances at the operating site.

⚠ **CAUTION:** For disinfection, if required, it is recommended to use water-based disinfectant solutions, with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated as they can discolour and/or damage plastic materials. This also applies to chemicals such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.

⚠ **CAUTION: Cleaning instruments.** Use appropriate cleaning instruments as specified in this manual. Never use brushes, metallic sponges, abrasive instruments or sharp objects to remove residues during manual cleaning. These may ruin the product surface and permanently and definitively damage the instruments. Only use brushes with soft, nylon bristles.

⚠ **CAUTION:** Use either tap water or demineralised water, as indicated by this manual.

⚠ **CAUTION:** Dispose of the enzymatic detergent correctly once used. Do not recycle.

3.1.3.3 Immersion in Enzymatic Detergent

This procedure applies to the following accessories:

- Handpiece;
- Torque wrench.
- Metallic front cone;

Prepare an enzymatic detergent solution
c) according to the manufacturer's instructions;


- c) Process validated by independent body with enzymatic detergent Enzymec (Mectron).

Place the handpiece, metallic front cone and torque wrench in a container. Add sufficient enzymatic detergent solution to completely cover the accessories;


⚠ CAUTION: Do not immerse the handpiece if warm.

Leave the handpiece, metallic front cone and torque wrench soaking for 10 minutes at 40 °C ±2; This will reduce the organic residues on the instruments;


1

2

3

Using a clean soft bristled, nylon brush under running water, clean thoroughly but without damaging the surfaces of the handpiece, front cone and torque wrench. Repeat until all visible traces of soil have been removed.

Also rinse the handpiece cord;

⚠ WARNING: The following parts must be brushed carefully:

Handpiece

- Threaded part onto which the inserts are screwed;
- Threaded part onto which the metallic front cone is screwed;
- Parts adjacent to threaded parts;

Metallic front cone

- The internal threaded part;

Torque wrench

- External metallic ring;
- Internal channels and cavities.




NOTE: Rinse the cavities, grooves and fissures and critical parts thoroughly.

⚠ CAUTION: To prevent the accumulation of enzymatic detergent residues, rinse the instrument surfaces thoroughly.


HANDPIECE:

Use a single-use syringe (min. 50 ml) to inject the enzymatic detergent solution into the cavities and internal channels of the handpiece three times in order to effectively remove the residues from the internal channel surfaces;

4


5

Rinse the internal channels, injecting cold distilled/demineralised water at a pressure of 3.8 bars for 10 seconds, to eliminate all enzymatic detergent residues.

NOTE: Rinse the channel within the handpiece by injecting the liquid from the back, taking care to use a waste water container to avoid sprays. Leave the water to drain through the internal channel. Repeat the procedure until the rinsing water is free from residues.

6

EN

3.1.3.4 Immersion in Neutral Detergent


This procedure applies to the following accessories:

- Handpiece;
- Metallic front cone.

Prepare a neutral pH cleaner solution^{d)} according to the manufacturer's instructions;


- d) Process validated by independent body with enzymatic detergent Enzymec (Mectron).

7

8

Place the handpiece and metallic front cone in a container. Add sufficient neutral pH cleaning solution into the container to completely cover the handpiece and front cone;

Use a clean soft bristled, nylon brush under running water to clean thoroughly but without damaging the surfaces of the handpiece and cone. Do this to remove all visible traces of contamination. Also rinse the handpiece cord;

⚠ WARNING: The following parts must be brushed carefully:

Handpiece

- Threaded part onto which the inserts are screwed;
- Threaded part onto which the metallic front cone is screwed;
- Parts adjacent to threaded parts



Metallic front cone

- The internal threaded part.


Rinse the internal channels, injecting cold distilled/demineralised water at a pressure of 3.8 bars, to eliminate all detergent residues.

NOTE: Rinse the channel within the handpiece by injecting the liquid from the back, taking care to use a waste water container to avoid sprays. Leave the water to drain through the internal channel. Repeat the procedure until the rinsing water is free from residues.

9

10

3.1.3.5 Ultrasonic Bath

This procedure applies to the following accessories:

- Torque wrench.

Prepare an enzymatic detergent solution^{e)} according to the manufacturer's instructions and pour it into an ultrasonic tank;

- e) Process validated by independent body with enzymatic detergent Enzymec (Mectron).


Place the torque wrench into a suitable tray (stainless steel washing machine tray with partitions).
Place the tray in the ultrasonic tank for 10 minutes at 40 °C;

Remove the torque wrench from the ultrasonic cleaner and rinse it with cold demineralised water;


NOTE: Rinse the cavities, grooves and crevices and all difficult-to-reach areas thoroughly.

Brush the torque wrench surface gently with a clean soft bristled, nylon brush.


11


12

13

14

3.1.4 Automatic Cleaning of Accessories

3.1.4.1 Washer - Disinfector

Place the handpiece, metallic front cone and torque wrench in a metallic tray. Arrange the material in a washer/disinfector into which the detergent solution and/or neutralising solution have been previously loaded, in accordance with the manufacturer's instructions;


NOTE: Place the instruments in the washer/disinfector so that dead zones do not arise and the water can properly drain. Also, make sure that the handpiece and the metal front cone are properly held in place in the washing tray and cannot move during the washing process, as shocks could damage them.

⚠ WARNING: Avoid overloading the washer/disinfector as this may compromise cleaning effectiveness.


⚠ CAUTION: Do not kink the handpiece cord; this may damage it.

Connect the specific adapter (supplied as optional) to the back of the handpiece and then to the appropriate water-jet connections of the washer/disinfector machine.

1

2


Sequence and parameters applicable to the cycle ^{f)}:

- 1 min, Pre-wash with cold water;
- 5 min, Washing with alkaline detergent at 55 °C $\pm$ 2 °C;
- 1 min, Neutralisation with suitable solution (1/3 cold e 2/3 warm);
- 1 min, Rinsing with tap water (1/3 cold e 2/3 warm);
- 5 min, Thermodisinfection at 93 °C with demineralized water.

The automatic disinfection was not experimentally tested. According to ISO 15883-1, the thermodisinfection at a temperature of 90 °C for 5 minutes determines an A0 value of 3000.

- f) Validated in compliance with standard ISO 15883-1 with Miele Washer - Disinfector, Program DES-VAR-TD; Alkaline detergent: Neodisher FA (0.2% v/v cleaning solution) and Neodisher Z (0.1% v/v neutralising solution).

3

EN


3.2 Cleaning Check

All instruments must be inspected before they are prepared for sterilisation. A visual inspection must be performed by means of a magnifier under a suitable light source.

All instrument surfaces must be controlled to reveal the presence of soils, corrosion and wear. Special attention must be paid to the control of critical points such as lumens, threaded parts, etc.

If soils should be seen, repeat the chosen cleaning cycle.


1

3.3 Drying and Lubrication of Accessories


Dry the handpiece, metallic front cone and torque wrench with a clean lint-free cloth;

1

Use filtered medical compressed air to effectively dry the handpiece, metallic front cone and torque wrench;

2

⚠ WARNING: Dry the following thoroughly:

Handpiece:


- Internal channel;

Torque wrench:


- Cavities, grooves and fissures and critical parts.

Lubricate the wrench using a medical lubricant for surgical instruments. Spray the lubricant directly on the internal peripheral contact surfaces of the torque wrench. After the lubrication remove the excess oil with a cloth.

3

⚠ CAUTION: It is important to lubricate the wrench and gearing after each cleaning cycle, only using a medical lubricant for surgical instruments in order to avoid any malfunction.

3.4 Sterilisation

3.4.1 Preparation

The following products require sterilisation:

- Handpiece with cord;
- Metallic front cone;
- Torque wrench.

The packaging must ensure the sterility of the instruments until opened for use at the sterile field.

The instruments must be packed with a double medical grade steam sterilization pouch, in compliance with the ISO 11607-1:2006 requirements.

⚠ CAUTION: Ensure that the inner pouch is large enough to contain the instruments without stressing the seals or tearing the packaging.

NOTE: The handpiece must be packed with the handpiece protection cap open.

NOTE: Each instrument must be packaged individually.

3.4.2 Materials Required

- Single-use pouches for steam sterilisation;
- Steam steriliser.

⚠ CAUTION: Sterilise only using a steam autoclave. Do not use any other sterilisation method as this may be incompatible with the materials used. Do not use the following sterilisation methods:

- Ethylene oxide sterilisation;
- Hot air sterilisation;
- Sterilisation with a flash cycle;
- STERRAD sterilisation;
- STERIS sterilisation;
- Gravity displacement autoclave;
- Other comparable sterilisation systems.

⚠ CAUTION: During the sterilisation cycle, do not exceed the allowed load of the steam sterilizer.

⚠ WARNING: Do not sterilise the instruments using:

- Hydrogen peroxide;
- Peracetic acid system;
- Sterilisation with formaldehyde vapours;
- Sterilisation with glutaraldehyde vapours;
- Other comparable sterilisation systems.

The sterilisation process must be carried out in a steam autoclave with a pre-vacuum cycle. Mectron S.p.A. guarantees SAL 10⁻⁶, by setting the parameters according to only one of the two validated procedures listed below.

Procedure A:

- **Cycle type:** 3 times pre-vacuum (min. pressure 60 mBar).
- **Minimum Sterilisation temperature:** 132 °C (range 0 °C ÷ +3 °C).
- **Sterilisation time:** 4 mins.
- **Drying time:** 10 mins.


or Procedure B

- **Cycle type:** 3 times pre-vacuum (min. pressure 60 mBar).
- **Minimum Sterilisation temperature:** 134°C (range -1 °C ÷ +1 °C).
- **Sterilisation time:** values equal to or between 3 mins and 18 mins.
- **Drying time:** 10 mins.

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

⚠ CAUTION: Do not sterilise the device console or pedal.

1



Sterile instruments must be stored in a clean, dry place at a storage temperature of between 5 °C and 40 °C.

! CAUTION:


- The handpiece and its cord cannot be separated.
- Do not sterilise the handpiece with the insert screwed into it.
- After sterilisation, the handpiece can only be used when it has cooled down to room temperature. The cooling process must not be accelerated.
- Re-screw the front cone to the handpiece before use.

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© Mectron S.p.A. 2021. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.