USE AND MAINTENANCE MANUAL

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EN

IT



DE

FR

ES

PIEZOSURGERY[®] flex

PIEZOSURGERY[®] flex

Copyright

© Mectron S.p.A. 2021. All rights reserved. No part of this document can be reproduced in any form without the written consent of the copyright owner.

EN

SUMMARY

1	Introduction	1
1.1	Intended Use	1
1.2	Description of the Device	2
1.2.1	Patient Group Directions	2
1.2.2	Patient Selection Criteria	2
1.2.3	Indications for Use	2
1.2.4	Users	3
1.3	Disclaimer	3
1.4	Safety Precautions	4
1.5	Symbols	6
2	Identification Data	8
2.1	Device Identification Label	8
2.2	Handpiece Identification Data	8
2.3	Inserts Identification Data	9
3	Delivery	10
3.1	List of the Components of the PIEZOSURGERY® flex	10
4	Installation	11
4.1	Safety Requirements in the Installation Phase	11
4.2	Connecting the Accessories	12
5	Use	16
5.1	Switching the Device On and Off	16
5.2	Description of the Keyboard	17
5.3	Foot Pedal Button	18
5.4	Safety Requirements Before and During Use	19
5.5	Instructions for Use	22
5.6	Important Information on Inserts	24
6	Maintenance	25
7	Disposal Method and Precautions	25
8	Technical Specifications	26
8.1	Electromagnetic Compatibility IEC/EN 60601-1-2	27
8.1.1	Guidance and Manufacturer's Declaration - Electromagnetic Emissions	27
8.1.2	Accessible Parts of the Casing	28
8.1.3	Guidance and Manufacturer's Declaration - Electromagnetic Immunity	29
8.1.3.1	Power Connection A.C. Input	29
8.1.3.2	Points of Contact with the Patient	31
8.1.3.3	Parts Accessible to the Input / Output Signals	32
8.1.4	Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device	33
9	Troubleshooting	34
9.1	Diagnostic System and Symbols on Keyboard	34
9.2	Troubleshooting Quick Guide	36
9.3	Fuses Replacement	38
9.4	Sending to an Authorized Mectron Service Centre	39
10	Warranty	40

PAGE INTENTIONALLY LEFT BLANK

PAGE INTENTIONALLY LEFT BLANK

1 INTRODUCTION

Carefully read this manual before proceeding with installation operations, use, maintenance or other intervention operations on the device. Always keep this manual at hand.

Important: To avoid damages to persons or things, read with particular attention all the chapters "Safety Requirements" present in the manual.

According to the degree of severity the safety requirements are classified with the following indications:

 **WARNING:** (Always referred to damage to persons)

 **CAUTION:** (Referring to possible damages to things)

The aim of this manual is to make the operators aware of the safety regulations, installation

procedures, instructions for proper use and maintenance of the device and its accessories. Use of this manual for aims other than those strictly linked to the installation, use and maintenance of the device is prohibited.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is engaged in continuously updating its products with possible changes to the components of the apparatus.

In case you encounter discrepancies between the descriptions found in this manual and the equipment in your possession you can:

- check for any available updates in the section MANUALS of MECTRON website¹;
- ask clarifications to Your Dealer;
- contact MECTRON After Sales Service.

1.1 Intended Use

PIEZOSURGERY® flex is a piezoelectric surgical system intended for use in bone cutting, equipped with special insert tips to perform osteotomy, osteoplasty and drilling in a variety of surgical procedures:

- Otolaryngology;
- Oral/Maxillofacial Surgery;
- Hand and foot Surgery;
- Plastic/Reconstructive Surgery.

It may also be used with endoscopic assistance to perform the above listed procedures.

 **WARNING:** Use the device only for the intended use. Failure to meet these requirements may cause serious injuries to the patient, the operator, and damages/breakdowns of the device.

 **WARNING:** Carefully read this manual and follow its recommendations in order to avoid to compromise the patient and/or the user safety. Failure to meet these requirements may cause serious injuries to the patient and/or the operator.

 **WARNING: Qualified and specialized personnel.** The device must be used exclusively by specialized personnel such as Surgeon. To use the device no special training is requested. The use of the device does not cause side effects if it is used correctly. An improper use might cause tissues heating.

 **WARNING:** The device must be used in a hospital environment, such as an operating theatre.

 **WARNING: Risk of explosion.** The device cannot function in places where there is an atmosphere saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

1.2 Description of the Device

The user interface has been optimized with the PIEZOSURGERY® flex making all the functions readily available by integrating them in the touch keyboard.

The PIEZOSURGERY® flex is a device that uses ultrasonic piezoelectric technology to generate mechanical microvibrations of the inserts, to effectively cut mineralized tissues. This allows an efficient and safe cutting which preserves the integrity of the osteotomized surfaces.

The micrometric, ultrasonic vibrations of the inserts provide greater precision and a selective cutting action compared to

traditional methods such as drills or oscillating saws (which act with macrovibrations), therefore minimizing traumatic effect on soft tissues.

The cavitation effect of the irrigating solution helps to keep the operatory field blood-free. This provides an optimal intra-operative visual control thus increasing safety, even in areas that are anatomically most difficult to access.

This medical device can be used on any patient of any age, weight, height, gender and nationality.

1.2.1 Patient Group Directions

This medical device is designed to be used with the following patient population:

- Infants;
- Children;
- Adolescents;
- Adults;
- Elderly.

This medical device can be used on any patient of any age, weight, height, gender and nationality.

1.2.2 Patient Selection Criteria

The use of the device is not recommended in the following cases:

1. Patients with active implantable medical devices (for example: pacemakers, hearing aids and/or other electromagnetic prostheses) without the prior authorization of their doctor;
2. Women who are pregnant or breastfeeding, due to the restrictions associated with the possible use of medical solutions such as anaesthetics
3. Patients with allergies;
4. Patients with pathologies or clinical conditions for which it is not recommended to perform surgery

or for which it may represent a contraindication according to the treating physician. These conditions may include, but are not limited to: heart disease, diabetes, cirrhosis, HIV infection, pregnancy or breastfeeding, radiation therapy, chemotherapy, immuno-suppressive therapy, allergies and psychiatric disorders;

5. Patients with unsuitable treatment sites.

All models of "Piezosurgery" bone surgery devices are intended for professional use only. Therefore, the user is the only person able to decide if and how to treat their patients.

1.2.3 Indications for Use

The use of the device is suitable for all the appropriate patients (see chapter above) for whom a bone surgery treatment is prescribed, by the treating physician, within the intended use of the device (see *Chapter 1.1 on page 1*).

1.2.4 Users

The device must be used only by specialised and properly trained personnel, such as a surgeon physician, able-bodied adults of any weight, age, height, gender and nationality.

1.3 Disclaimer

The manufacturer Mectron disclaim any liability, expressed or implied, and shall have no responsibility for any direct, indirect or other damages and personal injury arising out in connection with any improper practice in the use of the device and its accessories.

The manufacturer Mectron shall be under no liability, expressed or implied, with respect to any damages (personal injury and/or damage to property) which might arise or be caused, whether by the customer or by any of the users of the product and its accessories, as result of:

1. Procedures different than those specified in the intended use of the product;
2. The environmental conditions for the preservation and storage of the device are not compliant with the precautions indicated in the Chapter 8 on page 26;
3. The device is not used in compliance with all the instructions and precautions described in this manual;
4. The electrical system of the relevant operating room is not compliant with the applicable regulations and with electrical safety requirements;
5. The assembly operations, extensions, adjustments, updates, and repairs on the device are performed by personnel not authorized by Mectron;
6. Improper use, mistreatments, and/or incorrect interventions;
7. Any and all attempts to tamper with or modify the device, under any circumstance;
8. Use of non-original Mectron inserts that damage the threading of the handpiece, thus compromising correct operation and causing risk of harm to the patient;
9. Use of non-original Mectron inserts, even if they are used in accordance to designed and tested settings of Mectron original inserts. The correct use of the settings is guaranteed only with original Mectron inserts;
10. Lack of stock materials (handpiece, inserts, wrenches) to be used in the event of device stop due to fault or of inconveniences.

1.4 Safety Precautions

WARNING: Risk of explosion.

The device cannot operate in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

WARNING: Contraindications.

Interference with other equipment. The PIEZOSURGERY® flex device complies with the standard IEC 60601-1-2. However, it may interfere with other devices in its vicinity. Install PIEZOSURGERY® flex at safety distance from life-support systems. If adjacent or stacked use of the device is necessary, normal operation of the equipment and PIEZOSURGERY® flex, in the configuration in which they will be used, should be verified prior their use.

WARNING: Contraindications.

Interference by other equipment. An electrosurgical knife of other electrosurgical units placed in the vicinity of the PIEZOSURGERY® flex device may interfere with its correct functioning.

 **CAUTION:** In the case the hospital must subject the electro-medical equipment and systems to periodical inspections in order to adhere to imposed requirements in its own structure, the test procedures that must be applied to electro-medical equipment and systems to evaluate safety must be performed in line with norm EN 62353 'Electro-medical devices - Periodical inspections and tests to be performed after repair interventions on electro-medical devices'. The interval for periodic checks, in the intended operating conditions and described in this "Use and Maintenance" manual, is one year.

WARNING: Checking device status before the treatment.

Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. If anything unusual is noted during operation, do not perform the treatment. If the problem concerns the device, contact an Authorized Mectron Service Center.

 **CAUTION:** The electrical system of the relevant operating room, where the device is installed and used, must comply with the applicable regulations and with electrical safety requirements.

 **CAUTION:** To avoid risk of electric shock, this device must be grounded.

 **WARNING:** Do not operate the foot pedal of the PIEZOSURGERY® flex when the peristaltic pump cover is open. Moving parts could injure the operator.

 **WARNING: Personal injury.** Take care to ensure that cables do not hinder the free movement of personnel.

 **CAUTION:** US Federal Law (of the United States of America) restricts this device to sale by or on the order of a physician.

 **WARNING:** Prior to any use, system components should be inspected for damage. Do not use if damage is apparent.

 **WARNING:** Do not activate the device if the handpiece is faulty, damaged or broken. Replace the handpiece immediately.

 **WARNING:** Only use original Mectron accessories and spare parts.

 **CAUTION:** No modification of this equipment is allowed.

 **WARNING: Contraindications.** Do not use PIEZOSURGERY® flex on patients with heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ **CAUTION: Contraindications.** Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ **WARNING: Infection control.** To ensure maximum safety of the patient and the operator, verify the reusable parts and accessories has been cleaned and sterilized prior to use, following the instructions of the Cleaning and Sterilization Manual.

⚠ **WARNING: Cleaning and sterilisation of new or repaired tools.** All new and repaired reusable device accessories are supplied non-sterile. Before first use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilisation Manual.

⚠ **CAUTION: Contraindications.** When the reusable items (the handpiece, the torque wrench, and any other accessory that can be sterilized) have been autoclave sterilized, wait for them to cool down to room temperature prior to usage. The cooling process must not be accelerated.

⚠ **WARNING: Breakage and wear-out of the inserts.** In rare cases, high frequency oscillations and wear-out may lead to breakage of an insert. Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used. Should an insert fracture during use, extreme care must be exercised to ensure that all the fragments of the insert are retrieved and removed from the surgical site and, at the same time, that an effective suction is applied. During surgery, frequently check that the insert is intact, especially in its apical part. During surgery, avoid prolonged contact with retractors or other metal objects being used. During surgery, do not apply excessive pressure on the inserts.

⚠ **WARNING:** If an adverse event and/or serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the manufacturer indicated on the product label.

1.5 Symbols

Symbol	Description	Symbol	Description
	Device compliant with Regulation (EU) 2017/745. Notified body: IMQ S.p.A.		Nemko Mark UL CSA compliance
	Medical Device		Caution, read the Instructions For Use
	Operating Instructions		Manufacturer
	Date of manufacture		Serial Number
	Lot Number		Product Number
	Single-Use		Do not resterilize
	Use-by date		Sterilized with Ethylene Oxide (EO)
	Non-sterile		Can be sterilized in autoclave up to a maximum temperature of 135° C
	Type "B" applied part as per EN 60601-1 recommendation		Equipotentiality
	Alternating Current		Connection of the foot pedal
	Activation switch "on"		Activation switch "off"
	Functional earth		Caution, risk of electric shock
	Biohazard		The device and its accessories must not be disposed of or treated as solid urban wastes

Symbol	Description	Symbol	Description
	General warning ^{a)}		Do not touch moving parts ^{b)}
	Temperature limitation for transport and storage		Humidity limitation for transport and storage
	Atmospheric pressure limitation for transport and storage		Do not use if package is damaged
IPX8	International Protection Code of the mechanical casing. Liquid ingress protection.	IP20	International Protection Code of the mechanical casing.
QTY.1	Quantity of items in the package: 1	Rx Only	For US market, only. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed physician

Table 1 – Symbols


- a) The symbol is represented by a yellow warning triangle and a black graphical symbol.
- b) The symbol is represented by a red circle-with-slash and a black graphical symbol.

2 IDENTIFICATION DATA

A precise description of the model and device serial number will facilitate the After-Sales Service to answer the inquiries quickly and effectively.
Always provide this information when contacting an Authorised Mectron Service Centre.


2.1 Device Identification Label

Every device has an identification label that reports the main technical characteristics and the serial number. The identification label is located under the device. The complete technical specifications are reported in *Chapter 8* on page 26.

A separate label shows further symbols and characteristics of the device. This identification label is placed under the device.

NOTE: The complete list of symbols and their description are shown in *Chapter 1.5* on page 6.

2.2 Handpiece Identification Data

Each handpiece is laser-marked with (refer to *Figure 1* at page 9):

- its serial number (SN);
- its product number (REF);
- its date of manufacturing;
- the device name (PIEZOSURGERY® medical);
- the Mectron logo;
- HIBC 128 Data Matrix;
- CE, MD, Applied Part and Sterilizable symbols.

NOTE: The complete list of symbols and their description are shown in *Chapter 1.5* on page 6.



Figure 1 – Handpiece Identification Data

2.3 Inserts Identification Data

The following data are laser-marked on each insert:

- the name of the insert (Ref. 1);
- the Mectron logo (Ref. 2);
- the lot number to which the insert belongs (Ref. 3);
- the symbols “Single-use” and “CAUTION: read the instructions for use” (Ref. 4)

3 DELIVERY

3.1 List of the Components of the PIEZOSURGERY® flex

Refer to the front endpaper.

PIEZOSURGERY® flex consists of:

- A. Device core system;
- B. Peristaltic pump;
- C. Stand for Solution Bag;
- D. Use and Maintenance Manual, Cleaning and Sterilization Manual;
- E. Power Cord;
- F. Foot Pedal with bracket, cable and connector;
- G. Case.

Available accessories:

- H. PIEZOSURGERY® medical handpiece complete with cord and the connector protective cap (produced by Mectron);

⚠ CAUTION: The handpiece and the cord cannot be separated.

- I. PIEZOSURGERY® medical single-use inserts (produced by Mectron);
- J. PIEZOSURGERY® medical single use irrigation kit (distributed by Mectron);
- K. PIEZOSURGERY® medical torque wrench (produced by Mectron).

PIEZOSURGERY® flex consists of accessories that can be ordered separately.

The packaging of the device cannot undergo strong impacts as contains electronic components, therefore the transport and the storage must be carried out with particular care. All material shipped by Mectron is controlled at the time of dispatch.

The device is shipped appropriately protected and packaged.

Upon receipt of the device, check for any possible damage caused during transport and in case any damage and/or defects is found, complain to the transporter.

Keep the packaging in case there is a necessity to send any item to a MECTRON Authorized Service Centre and to store the device during long periods of inactivity.

⚠ WARNING: Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

4 INSTALLATION

To guarantee proper device functionality, it must be installed by an authorised Mectron technician. The device must be installed in a suitable and convenient place for its use. The technician will need to carry out the following 6 steps:

1. Unpack the device;
2. Explain to the user what precautions need to be taken to ensure correct installation;
3. Explain to the user all the possible device configurations;
4. Explain how to clean, sterilise and maintain the system;
5. Fill in the installation and hospital personnel training forms;
6. Send Mectron the completed forms to guarantee traceability of the device and the activation of the warranty.

4.1 Safety Requirements in the Installation Phase

⚠ WARNING: Contraindications. Interference with other equipment. The PIEZOSURGERY® flex device complies with the standard IEC 60601-1-2. However, it may interfere with other devices in its vicinity. Install PIEZOSURGERY® flex at safety distance from life-support systems. If adjacent or stacked use of the device is necessary, normal operation of the equipment and PIEZOSURGERY® flex, in the configuration in which they will be used, should be verified prior their use.

⚠ WARNING: Contraindications. Interference by other equipment. An electrosurgical knife of other electrosurgical units placed in the vicinity of the PIEZOSURGERY® flex device may interfere with its correct functioning.

⚠ WARNING: Risk of explosion. The device cannot operate in environments where the atmosphere is saturated with flammable gases (anesthetic mixtures, oxygen, etc.).

⚠ CAUTION: The electrical system of the relevant operating room, where the device is installed and used, must comply with the applicable regulations and with electrical safety requirements.

⚠ CAUTION: To avoid risk of electric shock, this device must be grounded.

⚠ WARNING: Do not operate the foot pedal of the PIEZOSURGERY® flex when the peristaltic pump lid is open. Moving parts could injure the operator.

⚠ WARNING: Install the device in a place protected against impacts or accidental splashes of water or other liquids.

⚠ WARNING: Do not install the device on top of or close to sources of heat. When installing, make sure that adequate air circulates around the device. Leave sufficient room around it, with particular reference to the fan on the rear.

⚠ CAUTION: Do not expose the device to direct sunlight or sources of UV light.

⚠ CAUTION: The device is transportable but must be handled with care when moved. Position the footswitch on the floor in such a way that it can only be activated intentionally by the operator.

⚠ CAUTION: Before connecting the handpiece cord to the device, make sure that the electrical contacts are perfectly dry. If necessary, dry them with compressed air.

⚠ CAUTION: Do not allow the device console and the foot pedal to get wet. If liquid enters the console or footswitch, damage could occur.

⚠ CAUTION: No modification of this equipment is allowed.

⚠ **CAUTION:** Position the device in way so that the power plug is always, easily reachable, since this plug is considered as a disconnecting means.

⚠ **WARNING:** Before using the PIEZOSURGERY® Medical irrigation kit, check the integrity of the sterile packaging and inspect the product to exclude the presence of any damage. Do not use the irrigation kit if the packaging is open or damaged. The irrigation set loses sterility if the packaging is broken or damaged. If the packaging is damaged or unintentionally opened before the use, proceed with the disposal of the set. The product cannot be re-sterilized or reused.

4.2 Connecting the Accessories

Equipotential plug: The device is equipped with an additional equipotential plug located on the rear of the console.


This plug is in accordance with DIN 42801.

Insert the connector of the equipotential cord (optional) to the plug on the rear of the device's console.


The purpose of additional potential equalization is to reduce differences of potential which can occur during operation between the device's body and conductive parts of other objects within the medical environment;

Plug the power supply cord into the power socket on the back of the device and then into the mains power outlet;

1


2

NOTE: The foot pedal comes equipped with a bracket that allows it to be moved to the place most suitable for the operation, without the need to use your hands to move it.

The bracket can also be positioned horizontally if it is not used;

3


Connect the foot pedal to the back of the device in the socket marked with the symbol , by means of the plug of the pedal cable, until you hear a "click" sound.

In order to disconnect the foot pedal from the device grab the connector of the foot pedal, press the release flap and pull the connector back;


NOTE: The button on the back of the device, left hand side, under the symbol , can be used to bring the treatment to its end in case the foot pedal is not functioning.

4

Fit the drip stand for supporting the solution bag into the dedicated hole;

5


STEPS TO BE CARRIED OUT IN THE "STERILE AREA":

Open the sterile package of the handpiece and of the irrigation kit, remove the tube and the clips;

Connect the end of the irrigation tubing onto the dedicated handpiece irrigation nozzle;

Use the 6 clips provided in order to clip together the irrigation tubing and handpiece cable;


6

STEPS TO BE CARRIED OUT IN THE "NOT STERILE AREA":


- Open the pump lid completely;

7

- Insert the thicker and 15 cm long tubing segment of the sterile line into the pump;
- Close completely the irrigation pump lid;

8

⚠ WARNING: When the peristaltic pump lid is open, the foot pedal of the PIEZOSURGERY® flex **MUST NOT** be activated. The parts in movement may cause harm to the operator.


Hang the saline bag to its specific drip stand;

Remove the protective cap from the spike and connect it to the irrigation bag;

⚠ WARNING: The irrigation bag support rod must be used with bags of up to a maximum of 1000 ml.

⚠ WARNING: The PIEZOSURGERY® Medical irrigation kit is provided in sterile packaging. Inspect the packaging to check its integrity. Whether it results damaged, DO NOT use the product and proceed with the proper disposal.


9

EN

Open the air inlet on the tube before operating. Release the irrigation tube clamp, if this is closed;


10

With the marked dot facing upward, insert the connector of the handpiece cord into the handpiece socket on the console. Gently push the handpiece-cord connector until it bottoms out.

⚠ CAUTION: To avoid damaging the handpiece-cord, grasp only the connector when disconnecting the cord. Never pull the cord.

11

5 USE

5.1 Switching the Device On and Off

Switching the Device On

Facing the front of the device, press the power switch on the left of the device body to the "I" position, being careful not to press the foot pedal. The device does a self-test and 4 symbols appear on the device (ref. P inside cover). At the end of the self-test the symbols turn off one by one, the device sets on the default setting and is ready to be used.


Switching the Device Off

Facing the front of the device, press the power switch on the left of the device body to the "O" position, being careful not to press the foot pedal. The device turns off.

NOTE: whenever the device is started, the following default setting is programmed:

- "Power" 1;
- "irrigation" 1;
- "mode" 1.


⚠ CAUTION: Once turned off, It is necessary to wait 5 seconds before switching the device on again. The device can give an error message.

5.2 Description of the Keyboard

TOUCH KEYBOARD

The user can configure the device by simply touching on the touch keyboard. Depending on the selected setting, the electronic feedback system automatically adjusts the correct operating frequency.

POWER (Ref. L front endpaper)


The desired level of power is adjustable from 7 levels, from 1 to 7, by selecting the numbers on the touch keyboard in the "power" column.

IRRIGATION (Ref. M front endpaper)

The delivery rate of the peristaltic pump is adjustable from 6 levels, by selecting the numbers on the touch keyboard in the "irrigation" column:


From 1 to 6= the pump flow goes from 8 ml/min to approximately 65 ml/min.

MODE (Ref. N front endpaper)

Depending on the type of surgery, it is possible to choose one of the 3 options available from the "mode" list:


1. Dedicated to the most delicate surgeries and to the lift of the sinus membrane;
2. Dedicated to the cut and to the removal of the mineralized bone;
3. Dedicated to the cut and to the removal of very thick mineralized bone.

FILLING THE IRRIGATION TUBING KIT (Ref. O front endpaper)

The device is equipped with the "pump" key which allows to perform the PUMP function.


The PUMP function can be used at the beginning of the treatment, to fill the entire irrigation tubing up to the insert, so that the surgery can be started with the necessary irrigation (see Chapter 5.5 on page 22).

SYMBOLS (Ref. P front endpaper)


PIEZOSURGERY® flex is equipped with a diagnostic system able to detect operating anomalies. Icons are displayed on the touchpanel according to the detected operating anomalies.

To help the user identify the malfunctioning part, four symbols are foreseen which are described in Chapter 9.1 on page 34.

5.3 Foot Pedal Button

In case the foot pedal is not functioning, the button on the back of the device, left hand side, under the symbol , can be used to bring the treatment to its end.

⚠ CAUTION: In the event you need to use the foot pedal button, first disconnect it from the socket of the unit.

⚠ WARNING: Use the foot pedal button only upon user request. The users must train their staff on how and when they have to use this button.

⚠ WARNING: The foot pedal button must be used only in case of malfunction of the foot pedal supplied with the unit. This button allows to end the surgical treatment in the event of a foot pedal malfunction.

5.4 Safety Requirements Before and During Use

⚠ WARNING: Before starting the activity, always assure having a set of spare parts (handpiece, inserts, wrenches) to be used in the event of a failure or problem.

⚠ WARNING: Only use original Mectron accessories and spare parts.

⚠ WARNING: Use of non-original Mectron inserts: this can completely damage the threading of the handpiece, affecting the proper operation and causing injury to the patient.

⚠ CAUTION: Contraindications. Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ WARNING: Contraindications. Do not use the PIEZOSURGERY® flex on patients with heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ WARNING: Checking device status before the treatment. Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. If anything unusual is noted during operation, do not perform the treatment. If the problem concerns the device, contact an Authorized Mectron Service Center.

⚠ WARNING: Infections control. First use: All new and repaired accessories are supplied in NON STERILE conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilization Manual. Subsequent uses: After every treatment, clean and sterilize all the reusable parts and accessories, following the instructions provided in the Cleaning and Sterilization Manual.

⚠ CAUTION: Contraindications. When the reusable items (the handpiece, the torque wrench, and any other accessory that can be sterilized) have been autoclave sterilized, wait for them to cool down to room temperature prior to usage. The cooling process must not be accelerated.

⚠ CAUTION: The electrical contacts inside the cord connector must be dry. Upon completion of the sterilization cycle, before connecting the handpiece cord to the device, make sure that the electrical contacts are perfectly dry. If necessary, dry them with compressed air.

⚠ WARNING: To ensure that the handpiece cools down, do not operate it without the irrigation circuit correctly installed and filled. To fill the irrigation circuit, always use the PUMP function.

⚠ WARNING: Always check delivery of the irrigation before and during use. Make sure the fluid comes out from the insert. Do not operate the device if the water is not delivered or if the pump is defective.

⚠ WARNING: Check the level of the physiological solution inside the relative saline bag. Replace the irrigation bag with a new saline bag before it is empty.

⚠ WARNING: The PIEZOSURGERY® Medical irrigation kit is guaranteed for one treatment only. Remove and dispose of the sterile line according to local regulations related to the hospital waste.

⚠ WARNING: Before using PIEZOSURGERY® Medical Irrigation Kit inspect the sterile package and the product for any damage. Do not use the irrigation kit if package is opened or damaged. Irrigation kit loses sterility in case of broken or damaged packaging. In case of damaged packaging, dispose of the irrigation kit. Do not sterilized or re-use the kit.

⚠ WARNING: Verify that the irrigation clamp is in the 'Open' position before operation, and that it is closed before disconnecting the irrigation tubing from the bag at the end of the procedure.

⚠ CAUTION: Never force the connector into the console port as this may damage the connector and/or console. If the connector and the port do not join with reasonable ease, they probably don't match. Make sure that the marked dot on the handpiece-cord connector is facing upwards.

⚠ CAUTION: Check that the PIEZOSURGERY® Medical handpiece is correctly connected before using the system.

⚠ WARNING: Before every treatment, make sure that the insert appropriate for the treatment is securely attached to the handpiece. Use exclusively the PIEZOSURGERY® Medical torque wrench to securing the insert to the handpiece. Do not use any other tool such as pliers, pincers, etc.

⚠ WARNING: Before every treatment, make sure that the selected insert has been correctly screwed onto the handpiece. This happens when the PIEZOSURGERY® Medical torque wrench used to tighten the insert emits a mechanical "CLICK" sound.

⚠ CAUTION: Hold only the connector when connecting/disconnecting the foot pedal. Never actually pull on the cord itself.

⚠ CAUTION: Do not attempt to screw or twist the connector during insertion or removal: the connector could get damaged by twisting.

⚠ CAUTION: The foot pedal is specifically designed to be used only in connection with the PIEZOSURGERY® flex device. Only use an original foot pedal otherwise damages or malfunctions can happen.

⚠ CAUTION: Do not activate the handpiece while the insert is in contact with the part to be treated. Doing so will not allow the electronic control circuit of the console to recognise the best point of resonance of the insert required for efficient, optimum performance.

⚠ WARNING: The patient must not come into contact with the device body or the foot pedal.

⚠ WARNING: Single use - Before the surgery. Before using any PIEZOSURGERY® Medical insert inspect the sterile package and the product for any damage. Insert loses sterility in case of broken or damaged packaging. In case of damaged packaging, the insert **MUST** be disposed of. **DO NOT** reprocess or reuse the insert whose packaging is damaged". Before starting the surgery correctly tighten the insert onto the handpiece by means of the torque wrench.

⚠ WARNING: Do not change the insert while the handpiece is operating to prevent possible user injury.

⚠ WARNING: During the tightening and removing operations, the user must pay particular attention to avoiding injury from inserts with sharp points and cutting edges.

⚠ WARNING: Breakage and wear-out of the inserts. In rare cases, high frequency oscillations and wear-out may lead to breakage of an insert. Do not bend, reshape or re-sharpen an insert in any way. Bending or prying the insert tip may cause it to fracture. Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used. During surgery, excessive pressure applied on the insert tip may cause it to fracture. If the insert tip contacts with metal objects, it may break, leaving pieces in the surgical site. Should an insert fracture during use, extreme care must be exercised to ensure that all the fragments of the insert are retrieved and removed from the surgical site and, at the same time, that an effective suction is applied. During surgery, frequently check that the insert is intact, especially in its apical part. During surgery, avoid prolonged contact with retractors or other metal objects being used.

⚠ WARNING: The PIEZOSURGERY® Medical device is intended for bone cutting. However, prolonged contact and/or excessive force of the instrument tip on soft tissues should be avoided as this may cause thermal and/or blunt injury. Particular care should be exercised when sharp tip inserts are used. Prolonged mechanical action of a sharp insert may also result in the soft tissue being cut. In close proximity to soft tissues/nerves (e.g. the perinevrium of the peripheral nerve system, or dura mater of central nerve system), it is recommended to complete the cut with a blunt tip diamond coated insert, to minimize the potential for soft tissue damage.

⚠ WARNING: Carefully check that the PIEZOSURGERY® Medical handpiece is correctly working in all its parts before using it on the patient.

ⓘ CAUTION: Intermittent operation. Long operational periods could cause overheating of the handpiece body. Refer to the Chapter 8 on page 26 for the recommended duty cycle.

⚠ WARNING: Sterility. The insert tip is supplied sterile by prior exposure to ethylene oxide gas (ETO).

⚠ WARNING: Single-use. The PIEZOSURGERY® Medical insert is intended to be used on an individual patient during a single surgical procedure and then discarded. The insert is not intended to be reprocessed, re-sterilized and used again. Remove and discard the insert following local regulations for proper disposal of contaminated surgical materials.

⚠ WARNING: During surgery, do not perform any maintenance activity on the system.


⚠ WARNING: The device is classified for intermittent mode operation (as reported in *Chapter 8 on page 28*). The continuous use of the device for a prolonged time and in any case exceeding the declared time limits can cause overheating, in particular, of the handpiece. In case of overheating, avoid contact of the handpiece with the operator and the patient.

5.5 Instructions for Use

After having connected all the accessories as described in *Chapter 4.2* on page 12, proceed as follows:

Screw the chosen insert onto the Piezosurgery® Medical handpiece till it bottoms out;


1

Use the Mectron torque wrench to secure the insert;


To use the Mectron torque wrench correctly, operate as follows:

2

Fit the insert into the wrench as shown;

3

Firmly hold the central body of the handpiece;

⚠ CAUTION: Do not grasp the end and/or cord of the handpiece; only hold the central body. Do not rotate the handpiece; simply hold it firmly and rotate the wrench, only.

Turn the wrench clockwise until it clicks (the external body of the wrench turns with respect to the body of the handpiece, making "CLICK" sounds).


The insert is now properly tightened;

Fill the irrigation tubing by selecting the PUMP function on the touch keyboard. The irrigation circuit starts to fill up;

As soon as the peristaltic pump starts, the entire scale of values of the "irrigation" column lights up and during the liquid passage the value of the irrigation shifts from 6 to 1;

As soon as the irrigation comes out from the Piezosurgery® Medical handpiece, the cycle can be stopped by pressing PUMP again or, alternatively, by pressing the foot pedal. The Pump function is disabled and the keyboard is enabled again, and displays the last setting used;

4


5

6


7

On the keyboard, select the necessary power, irrigation and mode;

⚠ CAUTION: According to the selected insert, for the correct setting of the parameters consult the Chart annexed to this manual titled "Recommended settings for the inserts" on the Piezosurgery® flex or the illustrative leaflet of the Piezosurgery® Medical insert you've purchased.

8

5.6 Important Information on Inserts

⚠ WARNING:

- The insert must always be kept moving. If the insert is blocked, the bone may overheat. We recommend using a continuous movement to minimise contact between the tip and the tissue. We recommend using high levels of irrigation as the power level is increased.
- For best performances, apply moderate, constant pressure on the insert. Do not apply excessive pressure; let the ultrasonic vibrations operate.
- Do not activate the handpiece while the insert is in contact with the part to be treated. Doing so will not allow the electronic control circuit of the console to recognise the best point of resonance of the insert required for efficient, optimum performance.
- Use original Piezosurgery® Medical inserts only. Use of non-original inserts, in addition to voiding the warranty, damages the threading of the Piezosurgery® Medical handpiece, with the risk of no longer being able to screw the original inserts correctly during subsequent use. Moreover, the device settings are tested and guaranteed to operate correctly only when original Piezosurgery® Medical inserts are used.

- Do not reshape, bend or re-sharpen an insert in any way. Doing so may cause the insert to fracture.
- Deformed or otherwise damaged inserts must never be used.
- Always check that the threaded parts of the insert and of the handpiece are perfectly clean – see the Cleaning and Sterilization Manual.
- During surgery, excessive pressure applied on the insert tip may cause it to fracture, resulting in harm to the patient.
- For information on how to correctly use the inserts, consult the Chart annexed to this manual titled "Recommended settings for the inserts" on the Piezosurgery® flex or the illustrative leaflet of the Piezosurgery® Medical insert you've purchased.
- Before using the PIEZOSURGERY® flex, make sure you have prepared the surgical site by moving away the soft tissues, to avoid damaging them. While cutting the bone, an accidental contact of certain parts of the insert with the soft tissues may cause them small injuries. To minimize this risk, use specific protective instruments.

6 MAINTENANCE

If the device is not used for a long time, observe the following recommendations:

1. Disconnect the device from the mains;
2. In case of a long period of inactivity, place the device in its original packaging and in a safe area;
3. Before using the device again, clean and

sterilise the handpiece and the wrench according to the instructions provided in the Cleaning and Sterilisation Manual;

⚠ WARNING: Regularly check that the power supply cord is in good condition. If damaged, replace it with an original Mectron spare part.

7 DISPOSAL METHOD AND PRECAUTIONS

⚠ WARNING: Hospital waste.

Treat the following items as hospital waste:

- Inserts at the end of each surgery;
- Irrigation kit at the end of each surgery;
- Torque wrench: when worn out or broken.

Disposable materials and materials that imply a biological risk  must be disposed of in accordance with current local regulations governing hospital waste.

The PIEZOSURGERY® flex must be disposed of and treated as a waste requiring separate collection.

The customer is entitled to deliver the old device for disposal to the retailer supplying the new equipment. Instructions on correct disposal are available from Mectron.

Non-compliance with the previous points may produce a fine in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive.

8 TECHNICAL SPECIFICATIONS

Device compliant with Regulation (EU) 2017/745	Classe IIa
Classification as per IEC/EN 60601-1	I Applied part type B (insert) IP 20 (device) IP X8 (foot pedal model FS-06)
Essential performance	According to the standard IEC 80601-2-60 the device has no essential performance.
Device for intermittent operation	60sec. ON - 30sec. OFF with irrigation
Power supply voltage	100-240 V ~ 50/60 Hz
Max. power consumption	120 VA
Fuses	Type 5 x 20 mm, T 2AL, 250V
Operating frequency	Automatic scan From 24 KHz to 36 KHz
Power settings	Adjustable on the touch screen: 7 power levels, from 1 to 7
Modes	Adjustable on the touch screen: from 1 to 3
Peristaltic pump capacity	Adjustable on the touch screen: 6 flow levels, from 1 to 6 (from 8 ml/min to approximately 65ml/min)
APC circuit protections	No handpiece detected Cord interruption Insert not tightened correctly or broken
Operating conditions	from 10 °C to 35 °C Relative humidity from 30% to 75% Air pressure P: 800hPa/1060hPa
Transport and storage conditions	from -10 °C to 60 °C Relative humidity from 10% to 90% Air pressure P: 500hPa/1060hPa
Altitude	Less than or equal to 2000 meters
Weight and size	3,2Kg 300 x 250 x 95 mm (L x l x H) ^{a)}

Table 2 – Technical Specifications

a) W = width ; L = length ; H = height

8.1 Electromagnetic Compatibility IEC/EN 60601-1-2

⚠ WARNING: Contraindications. Interference with other devices. The PIEZOSURGERY® flex device complies with the standard IEC 60601-1-2. However, it may interfere with other devices in its vicinity. Install PIEZOSURGERY® flex at safety distance from life-support systems. If adjacent or stacked use of the device is necessary, normal operation of the equipment and PIEZOSURGERY® flex, in the configuration in which they will be used, should be verified prior their use.

⚠ WARNING: Portable and mobile radio communication devices may influence the correct functioning of the device.

⚠ WARNING: Contraindications. Interference from other devices. An electro-surgical knife or other electro-surgical units placed in the vicinity of the PIEZOSURGERY® flex device may interfere with its correct functioning.

⚠ WARNING: The device requires specific EMC precautions and must be installed and commissioned according to the EMC information provided in this chapter.

⚠ WARNING: The use of cables and accessories not supplied by Mectron might negatively affect the EMC performances.

8.1.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

PIEZOSURGERY® flex is intended for use in the electromagnetic environment specified below. The customer or user of PIEZOSURGERY® flex should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	PIEZOSURGERY® flex only uses RF energy for internal function. Therefore, its FR emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF Emissions CISPR 11	Class B	PIEZOSURGERY® flex is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Compliant	

8.1.2 Accessible Parts of the Casing

PIEZOSURGERY® flex is designed to operate in the electromagnetic environment specified below.

The customer or user of PIEZOSURGERY® flex should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30 %.
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80 % AM a 1 kHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.
RATED power frequency magnetic fields ^{d) e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz o 60 Hz	The magnetic fields at the mains frequency should have levels characteristic of a typical location in a commercial or hospital environment.

- a) The interface between the PATIENT physiological signal simulation, if used, and the PIEZOSURGERY® flex shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the PIEZOSURGERY® flex.
- b) PIEZOSURGERY® flex that intentionally receives RF electromagnetic energy for the purpose of its operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.
- c) Testing may be performed at other modulation

frequencies identified by the RISK MANAGEMENT PROCESS.

- d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.
- e) During the test, the PIEZOSURGERY® flex may be powered at any NOMINAL input voltage, but with the same frequency as the test signal.
- f) Before modulation is applied.
- g) This test level assumes a minimum distance between the PIEZOSURGERY® flex and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that PIEZOSURGERY® flex will be used closer than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.

8.1.3 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

8.1.3.1 Power Connection A.C. Input

PIEZOSURGERY® flex is designed to operate in the electromagnetic environment specified below.

The customer or user of PIEZOSURGERY® flex should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic Environment Guidance
Electrical fast transients / bursts ^{a) l) o)}	IEC 61000-4-4	±2 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Surges ^{a) b) j) o)} Line-to-line	IEC 61000-4-5	± 0.5 kV, ± 1 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Surges ^{a) b) j) k) o)} Line-to-ground	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields ^{c) d) o)}	IEC 61000-4-6	3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in the ISM bands between 0.15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 KHz ^{e)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.
Voltage dips ^{f) p) r)}	IEC 61000-4-11	0% UT; 0,5 cycle ^{g)} At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	The quality of the network voltage should be that of a typical commercial or hospital environment.
		0 % UT; 1 cycle ^e 70 % UT; 25/30 cycle ^{h)} Single phase: at 0°	
Voltage interruptions ^{f) i) o) r)}	IEC 61000-4-11	0% UT; 250/300 cycle ^{h)}	The quality of the network voltage should be that of a typical commercial or hospital environment.

a) The test may be performed at any one power input voltage within the PIEZOSURGERY® flex RATED voltage range. If the PIEZOSURGERY® flex is tested at one power input voltage, it is not necessary to re-test at additional voltages.

b) All PIEZOSURGERY® flex cables are attached during the test.

c) Calibration for current injection clamps shall be performed in a 150 Ω system.

- d) If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) ME EQUIPMENT and ME SYSTEMS with a d.c. power input intended for use with a.c.-to-d.c. converters shall be tested using a converter that meets the specifications of the MANUFACTURER of the ME EQUIPMENT or ME SYSTEMS. The IMMUNITY TEST LEVELS are applied to the a.c. power input of the converter.
- g) Applicable only to PIEZOSURGERY® flex connected to single-phase a.c. mains.
- h) E.g. 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). PIEZOSURGERY® flex with battery backup shall resume line power operation after the test. For ME EQUIPMENT and ME SYSTEMS with RATED input current not exceeding 16 A, all phases shall be interrupted simultaneously.
- j) ME EQUIPMENT and ME SYSTEMS that does not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and ± 1 kV line(s) to line(s).
- k) Not applicable to CLASS II PIEZOSURGERY® flex.
- l) Direct coupling shall be used.
- m) r.m.s. , before modulation is applied.
- n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- o) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.
- p) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase.
- q) At some phase angles, applying this test to ME EQUIPMENT with transformer mains power input might cause an overcurrent protection device to open. This can occur due to magnetic flux saturation of the transformer core after the voltage dip. If this occurs, the ME EQUIPMENT shall provide BASIC SAFETY during and after the test.
- r) For ME EQUIPMENT and ME SYSTEMS that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the minimum and maximum RATED input voltage. ME EQUIPMENT and ME SYSTEMS with a RATED input voltage range of less than 25 % of the highest RATED input voltage shall be tested at one RATED input voltage within the range.

8.1.3.2 Points of Contact with the Patient

PIEZOSURGERY® flex is designed to operate in the electromagnetic environment specified below.

The customer or user of PIEZOSURGERY® flex should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{c)}	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30 %.
Conducted disturbances induced by RF fields ^{a)}	IEC 61000-4-6	3 V ^{b)} 0.15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 KHz	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

a) The following apply:

- All PATIENT-COUPLED cables shall be tested, either individually or bundled
- PATIENT-COUPLED cables shall be tested using a current clamp unless a current clamp is not suitable. In cases where a current clamp is not suitable, an EM clamp shall be used.
- No intentional decoupling device shall be used between the injection point and the PATIENT COUPLING POINT in any case.
- Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- Tubes that are intentionally filled with conductive liquids and intended to be connected to a PATIENT shall be considered to be PATIENT-COUPLED cables.
- If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test

frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

- The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18, 17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29, 7 MHz and 50,0 MHz to 54,0 MHz.
- b) R.M.S., before modulation is applied.
- c) Discharges shall be applied with no connection to an artificial hand and no connection to PATIENT simulation. PATIENT simulation may be connected after the test as needed in order to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

8.1.3.3 Parts Accessible to the Input / Output Signals

PIEZOSURGERY® flex is designed to operate in the electromagnetic environment specified below.

The customer or user of PIEZOSURGERY® flex should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{e)}	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30%.
Electrical fast transients / bursts ^{b) f)}	IEC 61000-4-4	±1 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Surges Line-to-ground ^{a)}	IEC 61000-4-5	± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields ^{b) d) g)}	IEC 61000-4-6	3 V ^{h)} 0.15 MHz - 80 MHz 6 V ^{h)} in the ISM bands between 0.15 MHz and 80 MHz ⁱ⁾ 80 % AM a 1 KHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

- a) This test applies only to output lines intended to connect directly to outdoor cables.
- b) SIP/SOPS whose maximum cable length is less than 3 m in length are excluded.
- c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) Calibration for current injection clamps shall be performed in a 150 Ω system.
- e) Connectors shall be tested per 8.3.2 and Table 4 of IEC 61000-4-2:2008. For insulated connector shells, perform air discharge testing to the connector shell and the pins using the rounded tip finger of the ESD generator, with the exception that the only connector pins that are tested are those that can be contacted or touched, under conditions of INTENDED USE, by the standard test finger shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) Capacitive coupling shall be used.

- g) If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- h) R.M.S., before modulation as applied.
- i) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18, 17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29, 7 MHz and 50,0 MHz to 54,0 MHz.

8.1.4 Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

PIEZOSURGERY® flex is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The customer or the user of PIEZOSURGERY® flex can help prevent electromagnetic interference by ensuring a minimum distance between the mobile and portable RF (transmitters) communication devices and PIEZOSURGERY® flex, as recommended, in relation to the maximum output power of radiocommunications equipment.

Test Freq. (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Max power (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 Band LTE 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2750	Bluetooth WLAN 802.11 b/g/n RFID 2450 Band LTE 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5420	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the PIEZOSURGERY® flex may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

⚠ WARNING: Portable RF communication equipment (including peripheral devices as antenna cables and external antennas) must not be used closer than 30 cm (12 inches) to any part of the device PIEZOSURGERY® flex, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these devices.

EN

9 TROUBLESHOOTING

9.1 Diagnostic System and Symbols on Keyboard

PIEZOSURGERY® flex is equipped with a diagnostic system able to detect operating anomalies. Icons are displayed on the touchpanel according to the detected operating anomalies. By using the following chart, the user is guided toward the identification and possible solution of the malfunction detected.

Symbol on the Keyboard	Possible Cause	Solution
	Handpiece electric contacts/ lead wet.	Thoroughly dry the contacts with compressed air.
	Piezosurgery® Medical handpiece not connected to the device.	Connect the handpiece.
	Faulty Handpiece.	Replace the handpiece.
	Malfunction of the tuning circuit.	Contact an Authorized MECTRON Service Centre.
	Insert not properly tightened on the handpiece.	Unscrew the insert and screw it back correctly using the torque wrench (see <i>Chapter 5.5 on page 22</i>).
	Insert broken, worn or deformed.	Replace the insert.
	Handpiece electric contacts/ lead wet.	Thoroughly dry the contacts with compressed air.

Symbol on the Keyboard	Possible Cause	Solution
	Peristaltic pump malfunction.	Check that there are no impediments to pump rotation.
	Irrigation kit not positioned correctly inside the pump.	Correctly reposition the silicon tubing inside the pump (see <i>Chapter 4.2</i> on page 12).
	The device has been turned off and on again without waiting 5 seconds.	Turn device off and wait 5 seconds before turning it on again.
	Abnormalities in the power distribution network, excessive electrostatic discharges or internal anomalies.	Turn the device off and wait 5 seconds before turning it on again. If the message persists, contact an Authorized Mectron Service Center.
	Turn-on procedure incorrect: the device has been turned on with the foot pedal pressed.	Check that during system start-up the pedal is not pressed. If the problem persists, disconnect the pedal and contact an Authorized Mectron Service Center.

Tabella 3 – Diagnostic System and Symbols on Keyboard

NOTE: For any error message not listed in the table, contact a Mectron Service Centre.

9.2 Troubleshooting Quick Guide

Problem	Possible Cause	Solution
The device does not turn on after having brought the switch into position "I"	The power supply connector is not properly plugged into the socket on the back of the device.	Check that the power supply cord is firmly connected.
	The power supply cord is faulty.	Check that the power supply socket is functioning. Replace the power supply cord.
	The fuses are out-of-order.	Replace the fuses (refer to <i>Chapter 9.3 on page 38</i>)
The device is switched on but not working. The display does not show any error.	The foot pedal plug is not correctly inserted into the device socket.	Correctly insert the foot pedal plug into the socket on the back of the device.
	The foot pedal is not working properly.	During the surgery: disconnect the foot pedal from the console and use the button replacing the pedal (see <i>Chapter 5.3 on page 18</i>). At the end of the intervention contact an Authorized Mectron Service Center
The device is switched on but not working. One of the following symbols appears on the screen: 	See <i>Chapter 9.1 on page 34</i> for the possible cause, according to the symbol that has been displayed.	See <i>Chapter 9.1 on page 34</i> for the action to undertake, according to the symbol that has been displayed
When operating, the Piezosurgery® Medical handpiece emits a slight whistle.	The insert is not properly tightened into the handpiece.	Unscrew the insert and screw it back correctly by using the Mectron torque wrench. (Refer to <i>Chapter 5.5 on page 22</i>)
	The irrigation circuit has not been completely filled.	Fill the irrigation circuit by using the "pump" function (Refer to <i>Chapter 5.5 on page 22</i>)

Problem	Possible Cause	Solution
During operation, no irrigation flows out from the insert.	The insert is clogged.	Unscrew the insert from the handpiece and unclog the insert water passage by blowing compressed air through it. If the problem persists, replace the insert with a new one.
	The handpiece is clogged.	Contact an Authorised Mectron Service Centre.
	The liquid irrigation bag is empty	Replace the empty bag with a new one full.
	The air inlet on the spike has not been opened.	Open the air inlet on the spike/roller.
	Irrigation kit not positioned correctly inside the pump.	Check that the Irrigation kit is properly connected.
The device works properly but the pump strains/ is stressed.	Excessive pressure of the impeller on the peristaltic pump tube.	Check that the tube inside the peristaltic pump is properly positioned (Refer to <i>Chapter 4.2 on page 12</i>)
The pump turns correctly, but when it stops, irrigation flows out from the handpiece.	The lid of the peristaltic pump is not closed correctly.	Check that the lid of the peristaltic pump is perfectly closed (Refer to <i>Chapter 4.2 on page 12</i>)
Poor performance	The insert is not correctly tightened on the handpiece	Unscrew the insert and screw it back correctly by using the Mectron torque wrench. (Refer to <i>Chapter 5.5 on page 22</i>)
	Broken, worn-out or deformed insert.	Replace the insert with a new one.

Tabella 4 – Troubleshooting Quick Guide


9.3 Fuses Replacement

⚠ WARNING: Switch the device off.
Always turn the device off by means of the main switch and disconnect it from the electrical power socket before proceeding.

EN

Apply leverage with a flat screwdriver, inserting its tip in the seat of the fuse-holder drawer located under the power supply socket;


1

Pull out the fuse-holder drawer;


⚠ WARNING: Replace the fuses, complying to the characteristics indicated in *Chapter 8 on page 26*.

2

Reinsert the drawer in its housing.

3

9.4 Sending to an Authorized Mectron Service Centre

In the event that technical assistance is needed on the device, contact one of the Authorized Mectron Service Centers or your Retailer. Do not attempt to repair or modify the device and its accessories.

Clean and sterilize all the parts that need to be sent to an Authorized Mectron Service Center, following the instructions provided in the Cleaning and Sterilization Manual supplied with the device.

Leave the sterilized parts in the packet, procedure that confirms the sterilization process has been performed.

The cleaning and sterilization requirements are in line with those in force concerning the safeguarding of health and safety on the workplace, as per DLgs 81/08 and its subsequent amendments, both laws of the Italian State.

In the case that the customer does not fulfil the requirements, Mectron reserves the right to charge him or her the cleaning and sterilization expenses, or to reject the goods received in unsuitable conditions, returning them to the customer, at his or her expense, for them to be correctly cleaned and sterilized. The device must be returned properly packaged, accompanied by all the accessories

and by a sheet bearing the:

- Data of the owner with telephone number;
- Product name;
- Serial number and/or lot number;
- Reason for goods returned / description of the malfunction;
- Photocopy of delivery note or purchase invoice of the device.

⚠ CAUTION: Packaging

Pack the device in its original packaging to prevent damages during transport.

Once the material has been received by the Authorized Mectron Service Center, qualified technical personnel will provide its evaluation to the given circumstance. The repair works will be performed only upon their prior acceptance by the end customer. For further details, contact the closest Authorized Mectron Service Center or your retailer.

Unauthorized repair works may damage the system and make the warranty void, and absolve Mectron of any responsibility for direct or indirect harm to persons or damages to things.

10 WARRANTY

Before being marketed, all MECTRON devices, are subjected to a thorough final check that verifies the full functionality.

MECTRON warrants its products, purchased new from a MECTRON dealer or importer, against defects in material and workmanship for:

- 2 YEARS (TWO) for the device, starting from the installation and commissioning date;
- 1 YEAR (ONE) for the handpiece together with its cord, starting from the installation and commissioning date.

The other accessories are not included in the warranty.

During the period of validity of the warranty, MECTRON undertakes to repair (or, at his free choice, replace) free of charge those parts of products that in their opinion prove being defective.

MECTRON disclaims any responsibility for direct or indirect damage to people or things, and the manufacturer's warranty and device approval is not valid in the following cases:

- the device is not used according to the intended use for which it is provided;
- the device is not used in accordance with all the instructions and requirements described in this manual;
- the electrical system of the places where the equipment is used do not comply with the laws in force and the related regulations;
- Assembly operations, extensions, re-adjustments, modifications, replacements and repairs are carried out by personnel not authorized by Mectron or in breach of what is provided in this manual also in regard to the origin of the authorised material;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 7 on page 34*.
- non-original Mectron inserts, accessories or spare parts have been used that may

affect the correct working of the device, and cause injury to the patient.

- Accidental breakage is caused during transport.
- Damage caused by incorrect use or negligence, or to connection to a power voltage other than the one foreseen.
- The cleaning and sterilisation procedures have not been properly performed.
- The warranty has expired

NOTE: The warranty period starts from the installation date, as written on the installation and test module.

The warranty is void in case the device has been tampered with, or repaired by not authorized personnel. See Chapter 9.4 on page 39 for the details relative to the procedure to send the device to an Authorized Mectron Service Center.

The expected service life of the device is 5 years, minimum.

The service life/duration does not define a limit of use; the service life of the device defines the period of time, subsequent to installation and/or commissioning, during which the original performances or, in any case, adequate for the intended use are guaranteed, without any degradation occurring such as to compromise its functionality and reliability.

The service life is a minimum qualitative target of the design, therefore, it is not excluded that single parts or components guarantee performances and reliability higher than those declared by the manufacturer.

The service life is intended in compliance with the maintenance plans provided for in this manual, it does not include components normally subjected to "wear" and it is independent of the warranty period: the service life period does not establish any implicit or explicit extension of the warranty period.

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© Mectron S.p.A. 2021. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.