Original Operating Instructions

N4 Impression

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dentalportal.info

Welcome

This document is designed and released for the following groups/individuals:

- End users
- Authorized resellers
- Authorized service technicians

Distribution or duplication of this document is only allowed by written consent of vhf camfacture AG. This includes reproduction by presentation and broadcasting.

This document is an original document.

It is subject to technical changes and errors.

Contents

1	About this document	4
1.1	Used symbols	4
1.2	Structure of safety notes	4
2	General safety instructions and operating regulations	4
2.1	General safety instructions.....	4
2.2	Operating regulations	6
2.2.1	Intended use	6
2.2.2	Controlling the machine through software	6
2.2.3	Maintenance and cleaning	6
2.2.4	Spindle	6
2.2.5	Unattended operation	7
2.2.6	Transportation and storage	7
3	Machine description	8
3.1	Connection panel on the backside.....	8
3.2	Working chamber door	8
3.3	Working chamber	9
3.4	Cooling liquid tank	9
3.5	Manufacturing software & CAM computer	10
3.6	Sound emission	10
3.7	Technical data	10
4	Installing the machine	11
4.1	Checking the scope of delivery.....	11
4.2	Choosing the installation site.....	11
4.3	Machine installation (scheme)	12
4.4	Installing the pneumatics	13
4.4.1	Overview service unit	13
4.4.2	Mounting the service unit on the machine	13
4.4.3	Installing the pneumatic hoses to the service unit.....	14
4.4.4	Setting the air pressure with the service unit	14
4.5	Filling the cooling liquid tank	14
4.6	Establishing the electric connection	15
4.7	Removing the transport lock.....	15
4.8	Connecting the CAM computer	16
4.9	Testing the machine	16
5	Running the machine	17
5.1	Switching on the machine	17
5.2	Overview job execution.....	17
5.3	Inserting / changing tools	18
5.4	Mounting and removing blocks	19
5.5	Using an optional abutment holder	20
5.5.1	Installing the PreFace® abutment holder	20
5.5.2	Demounting the PreFace® abutment holder	20
5.6	Mounting and removing PreFace® abutments	20
5.7	Executing jobs.....	21
5.8	Interrupting the job execution.....	21
5.9	Operation interruption or abortion of job execution	21
5.9.1	How to proceed in case of an operation interruption	22
5.9.2	How to proceed in case of a machine malfunction.....	22
5.9.3	How to proceed in case of a tool breakage.....	22
5.9.4	Procedure in case of a power failure.....	22
5.10	Emergency release of the working chamber door	22
6	Machine maintenance and cleaning	24
6.1	Definition of wear parts.....	24
6.2	Maintenance table	24
6.3	Inspection.....	25
6.4	Cleaning the working chamber and measuring key ..	25
6.5	Cleaning the collet chuck	26
6.6	Checking the hoses, cables and connections	27
6.7	Removing the cooling liquid tank completely from the machine.....	28
6.8	Re-inserting the cooling liquid tank	28
6.9	Exchanging the cooling liquid / cleaning the cooling liquid tank	29
6.10	Exchanging the pump.....	30
6.11	Cleaning the housing.....	30
6.12	Checking the service unit	31
6.12.1	Checking the water separator for condensate.....	31
6.12.2	Exchanging / cleaning the contaminated filter cartridge	31
6.13	Exchanging the main fuse	31
6.14	Calibrating the axes	32
6.15	Replacing the tool magazine inserts	34
7	Disposal	35
7.1	Disposing of the cooling liquid	35
7.2	Disposing of the machine	35
8	Troubleshooting	36
Index		39
EC Declaration of conformity		41

1 About this document

1.1 Used symbols

Instructions

- ★ Prerequisite
- General instruction
- M Specific manual action
- S Specific action in the software
- Result

Additional symbols

-  Cross reference
- List
-  Additional information
-  Important instructions without any danger for people or objects.
-  Information to make work more efficient.

1.2 Structure of safety notes

CAUTION

Type and source of hazard

Further explanations and consequences when ignoring the hazard.

- Instructions to avoid hazards.

The following signal words may occur in this document:

DANGER

DANGER indicates a hazardous situation which will result in death or serious injury.

WARNING

WARNING indicates a hazardous situation which can result in death or serious injury.

CAUTION

CAUTION indicates a hazardous situation which can result in minor or moderate injury.

NOTICE

NOTICE indicates a situation which can lead to material damage of the product or in the surrounding areas(s).

2 General safety instructions and operating regulations

2.1 General safety instructions

DANGER


Incorrect operation of the machine

-  ➤ **Before** installing, maintaining and operating the machine, read **all** documents provided with the machine.
- If it is unclear how to operate the machine in any way, do not use the machine and contact customer service.
- Ensure that every user has access to the operating instructions.
- Instruct every user on safe and proper machine handling.

Danger to life due to an electric shock

 If you come into contact with electrically charged parts, you can suffer from an electric shock. Water increases the risk significantly.

- Do not remove the housing of the machine.
- Only have qualified electricians work on any electric equipment.
- Ensure that an operational Residual Current Device/Ground Fault Circuit Interrupter is installed on the electric circuit of the machine.
- Run power cables so that they cannot be damaged by sharp edges.
- **Before** switching on the machine, check power cables for damage.
- **Before** unplugging the power cable, switch off the machine at the main power switch.

 ➤ In the following cases disconnect the machine from the electrical source immediately and prevent it from being restarted:

- When machine connections or electric cables are damaged
- When liquid is leaking
- **Before** you check or run electric cables
- Replace damaged cables with original manufacturer's spare parts.
- Never perform any troubleshooting while the machine is operating.
- Only have authorized service technicians repair the machine.
- Do not touch the machine and especially the cables with wet or damp hands.

- Check the environment of the machine and all accessible internal areas daily for leaked liquid and remove any liquids near or in the machine immediately.
- Never put any machines or devices which are powered by electricity under the machine.
- Do not put any objects onto the machine.

WARNING

Respiratory diseases when processing harmful materials

If you inhale harmful materials during their processing, your respiratory tracts may be damaged.

- Avoid materials which damage your health.

Danger to health if using the wrong cooling lubricant

Some cooling liquids can be a serious threat to your health and /or the environment.

- Only add the cooling lubricant Tec Liquid Pro to the cooling liquid. The mixing ratio is indicated on the bottle label.

Crushing hazard and risk of cutting injuries through moving machine parts

Through the moving axes and the rotating spindle you can suffer bruises and cuts.

- Only operate the machine when the working chamber door is closed and the safety interlock is activated.
- Store the key for the emergency release in a place where only authorized persons have access.
- Do not circumvent or deactivate safety devices of the machine.
- Check the machine regularly for damage, especially the safety devices.
- Have damaged safety devices repaired by an authorized service center unless stated otherwise in this document.
- Only use original manufacturer's equipment and original spare parts for the machine.
- Keep children and animals away from the machine.
- Do not remove the housing of the machine.

Service Mode: risk of cutting injuries and bruises as well as hazards through ejected debris

If you operate the machine in any "Service Mode" with the working chamber door open, the risk of injury is increased significantly!

- Operate the machine in "User" mode only unless you have been authorized by vhf to use other modes.
- Even if you are an authorized user, use the "Service Modes" only when necessary.

- While in any "Service Mode": Do not reach into the working chamber while the axes are moving or during machining.

- While in any "Service Mode": All persons within the reach of the machine must wear protective eye wear.

Hearing damage due to loud noise

If you are regularly exposed to loud machining noise, you may suffer from hearing loss and tinnitus,

- If loud noise cannot be avoided, wear ear protection during machining.

Risk of injuries through loose pneumatic components under air pressure when connections are open

Loose pneumatic components can move extremely fast and unpredictably and may cause injury.

- **Before** you run the pneumatic hoses, close the compressed air supply valve.
- **Before** you check the pneumatic hoses and pneumatic connections, set the air pressure to a minimum value.

- In case of defective machine connections and pneumatic hoses, disconnect the machine from the external compressed air supply and the electrical source.
- Contact customer service if connections are damaged or defective.

CAUTION

Trip, fall and slipping hazards

- Keep the working environment and installation site clean.

- Run cables in such a way that persons cannot trip over them.

Risk of cutting injuries and burns

If you touch tools or sharp edges on blanks or the machine, you may suffer from cuts. If you touch the hot spindle body or hot tools, you may suffer from burns.

- Wear gloves when you perform manual work at the machine or with blanks / tools.

Danger to health if handling the cooling lubricant incorrectly

- **Before** using the cooling lubricant, read the safety data sheet that was supplied with the cooling lubricant.
- When handling the cooling lubricant, **always** wear suitable protective clothing.
- **Always** store the cooling lubricant in the original container.

Reduced ability to act with insufficient lighting

In case of an insufficient lighting your judgement and/or your precision may be reduced.

- Make sure that the lighting in your working environment is sufficient.

Risk of injury in case of malfunctions caused by insufficient maintenance

If you do not maintain the machine as required, malfunctions may occur which can lead to injuries.

- Take note of the intervals and conditions mentioned in the maintenance table in the operating instructions. Carry out the respective maintenance tasks accordingly.

Health risk through constant malpositioning if your working environment is not sufficiently ergonomic

Over the long run, an improper or one-sided positioning can be a risk to your health.

- Set up an ergonomic working environment.
 - Ensure that the seat height and monitor position is ideal and the lighting is sufficient.
-

2.2 Operating regulations

If you violate the following regulations, you may lose your entitlement to benefits.

NOTICE

Machine damage if you violate these regulations

If you violate the following regulations, your machine may get damaged and / or cause damage in the surrounding area(s).

- Thoroughly follow all instructions and information in this section.
-

2.2.1 Intended use

The machine and the manufacturing software have been designed for the commercial creation of dental objects by specially trained persons. The objects require additional work before their use on patients.

- Only process materials that you can select in the manufacturing software.
- Only use the machine and manufacturing software in a commercial environment.
- Before creating jobs, verify if the objects being prepared may be utilized at the place of use according to local and / or national regulations or other authorized organizations or entities (e.g. professional associations, health authorities). In particular, verify if the material is approved for the machined object type and if the object type is designed in accordance with applicable regulations. Neither the manufacturing software nor the machine will inform you about possible regulatory infringements, but will execute jobs in accordance with the preferences and materials set by the user.
- Verify that each object type and each material in your jobs are authorized manufacturing materials. If mandated by local or national regulations, obtain relevant authorization from the responsible organizations or entities (e.g. professional associations, health authorities).
- Only manufacture objects which correspond to the object types that you can select in the manufacturing software. While you can import / manufacture any other objects as well, neither the manufacturing software nor the machine are designed for these other objects and should not be used in this way.
- Do not manufacture implants or parts of objects that are designed to have contact with implants. These parts include parts of two-part abutments which contain the connection geometry for the implant. Do not manipulate the connection geometry of prefabricated abutments (“prefab abutments”) and you must always check finished objects for accurate connection geometries (i.e. that connection geometries of finished jobs have not been damaged).

2.2.2 Controlling the machine through software

You control the machine through specially designed applications which are supplied with the machine.

- Always use the latest program version that officially supports your machine.
- Before installing or operating the machine, read the documentation for the applications.
- Ensure that your CAM computer meets all system requirements.

2.2.3 Maintenance and cleaning

Maintenance and cleaning is part of standard machine usage.

- Clean and maintain the machine as required. Only then can the machine reach a long service life.

2.2.4 Spindle

The spindle of your machine is a high-precision instrument.

- Do not use unbalanced tools at high rotational speeds. Such an imbalance puts a great strain on the spindle's ball bearings, which can cause the bearings to be damaged.
- When working in the working chamber, do not apply manual pressure against the spindle.

2.2.5 Unattended operation

If the machine runs unattended, the risk of material damage is increased.

- Only allow unattended operation of the machine to occur if the following conditions are met:
 - The national and local laws allow it.
 - The working chamber is completely clean.
 - Unauthorized users cannot access the machine.
 - The room in which the machine is located has an automatic fire detection system.

2.2.6 Transportation and storage

WARNING

Injuries caused by unsafe transportation

If you transport the machine unsafely, the machine may slip and cause injuries.

- Always transport **unpacked** machines individually and do not stack them.
- Ensure that only trained personnel transport the machine to and from the installation site.
- Ensure that the housing of the machine is completely closed.
- Always transport the machine in an upright position.
- Transport and position the machine with 2 people.
- Grab unpacked machines only at the left and right grips located at the bottom of the machine.

NOTICE

Short-circuit hazard when the machine is too cold

If the machine is transported from a cold environment into a warmer environment, a short circuit may occur caused by condensate.

- **Before** switching on the machine after transportation, ensure the following:
 - The ambient air has the allowed temperature.
 - The machine has the same temperature as the ambient air. This will take **at least** 48 hours.
 - The machine is completely dry.

- Ensure that the following conditions are met during the whole transportation and/or storage period:

- Ambient temperature (storage/transport): between -20 °C (-4 °F) and 60 °C (140 °F)
- Relative air moisture: max. 80 %, non-condensing

Preparing transportation or storage

Before transporting or storing your machine, the following preparations are necessary:

- M1. Remove all blanks from the working chamber.
- M2. Rinse the cooling liquid system.
- M3. Drain and clean the cooling liquid tank.
- M4. Clean the working chamber.
- M5. Install the transport lock. For this, see the corresponding steps on the supplement.
- M6. Ensure that the housing of the machine is completely closed.
- M7. Switch off the machine at the main power switch.
- M8. Disassemble the machine components by following the installation instructions in reverse order.
- M9. In case of overseas transport, take proper measures against corrosion.

Repackaging

To repack the machine after preparing its transportation or storage, the following steps are necessary:

- M1. If possible, use the original packaging. If the original packaging is not available, use a packaging of similar size and quality.
 -  [Original packaging is available from customer service.](#)
- M2. Pack the machine and its accessories securely.
- M3. Protect the packaging against slipping. If machines are properly packed and protected against slipping, they may be stacked.

3 Machine description


FIG. 1: FRONT VIEW OF THE N4 IMPRESSION

- [1] Connection panel on the backside (↗ Fig. 2)
- [2] Recessed grip for opening the working chamber door
- [3] Working chamber door
- [4] View window to the working chamber
- [5] Cooling liquid tank

With your N4 Impression you can process blanks in block shape of different materials to create high quality objects for the dental sector. The tools and blanks are constantly being cooled by a cooling liquid (water and special cooling lubricant) during machining.

You can find a list of the materials which you can process with the N4 Impression in the DentalCAM manufacturing software.

3.1 Connection panel on the backside


FIG. 2: CONNECTION PANEL ON THE BACKSIDE

- [1] Pneumatic connection
- [2] USB connection
- [3] Power connection 100 – 240 V AC, 50/60 Hz, including glass fuse T6,3A L250V
- [4] Main power switch

3.2 Working chamber door

The working chamber door locks the working chamber and protects the user from injuries during operation.

CAUTION

Crushing hazard when opening or closing the working chamber door

When you open or close the working chamber door, the moving working chamber door may crush your fingers. Objects on the machine may fall and cause injury or damage.

- When you open and close the working chamber door, use one hand and keep the other hand away from the machine.
- When you close the working chamber door, ensure your hands do not get caught between the door and the machine housing.
- Do not place objects on the machine.

NOTICE

Damaging of the machine when opening the working chamber door by force

When the machine is not supplied with power or when the axes are moving, the working chamber door remains locked. If you open the locked working chamber door by force, the machine may get damaged.

- Never open the working chamber door by force.
- **Before** you open the working chamber door, connect the machine to the electrical source and switch it on.
- Always close the cooling liquid tank first before opening the working chamber door.
- To open the working chamber door, put one of your hands into the recessed grip and then pull the working chamber door down.
- To close the working chamber door, put one of your hands under the open working chamber door and push it upwards.

3.3 Working chamber

The machine processes all blanks in the working chamber.


FIG. 3: THE WORKING CHAMBER OF THE N4 IMPRESSION

- [1] Tool magazines for up to 8 tools in total
- [2] Bellow
- [3] Blank holder and rotational axis
- [4] Spindle with collet chuck for picking up tools

Colors of the working chamber lighting

! If the working chamber lighting is insufficient, provide additional lighting.

Your N4 Impression illuminates the working chamber in different colors. The color will change depending on the state of the machine. You find the colors and respective machine statuses in the following table:

Colour	Status
Green	Machine is ready, working chamber door closed
White	Machine is ready, working chamber door open
Blue	Machine is executing a job
Red	A machine malfunction has occurred

3.4 Cooling liquid tank

NOTICE

Damaging of the machine when using a non-approved cooling lubricant or machining with no lubricant at all

The approved cooling lubricant protects the machine from damage caused by wet machining.

- Only use cooling lubricants that have been approved for your machine by customer service.
- Before wet machining, ensure that the proper amount of cooling lubricant has been added to the cooling liquid.

You fill the cooling liquid tank with the cooling liquid. The cooling liquid is an emulsion of water and the cooling lubricant **Tec Liquid Pro** (chapter 4.5, page 14).

After injection into the working chamber, the cooling liquid flows back into the tank. Due to this closed circuit, the loss of cooling liquid during processing is minimal.

A washable coarse filter mat is installed in the cooling liquid tank and protects the pump from dirt.


FIG. 4: THE COOLING LIQUID TANK OF THE N4 IMPRESSION

- [1] Coarse filter mat
- [2] Tank for the cooling liquid

- Only open the cooling liquid tank when the working chamber door is closed.
- Never open the cooling liquid tank while the machine is operating.
- Immediately wipe up leaking cooling liquid under the machine after opening the cooling liquid tank.
- Always close the cooling liquid tank completely.

3.5 Manufacturing software & CAM computer

To operate the machine, you must use a computer running Windows® (“CAM computer”) which is connected to the machine via the provided USB cable.

You can process manufacturing jobs with the CAM computer in 2 steps:

1. With DentalCAM you can create jobs with the objects that you want to machine. DentalCAM already contains all important parameters for the different materials.
2. You can execute the jobs with the CNC software DentalCNC. This application offers additional maintenance and control functions for the machine.

i You must use a different computer together with a scanner and suitable software to create the 3D models of the objects (not provided, sold by specialist dealers).

3.6 Sound emission

The actual sound emission of the machine varies heavily depending on the manufacturing material and the grinding conditions.

- If the machine is too loud, check the operating conditions. Ensure that the blank is fixed properly, check the condition of the tool and the material you use.
- If loud noise cannot be avoided, wear ear protection during machining.

Sound measurement

Measuring conditions:

- Processed material: VITABLOCS Mark II®
- State of the tool: new
- Distance to sound source: 1 m
- Measurement according to ISO 3746, engineering method 3

Operating condition	Maximum A-weighted sound pressure level
Processing	78 dB(A)
All other operating conditions (tool change, movement of the axes etc.)	<70 dB(A)

3.7 Technical data

Base system

- Dimensions (W/D/H):
 - Footprint: approx. 340 x 320 mm (13.4 x 12.6 in)
 - Housing fully closed: approx. 360 x 450 x 470 mm (14.2 x 17.7 x 18.5 in)
 - Housing fully open: approx. 360 x 690 x 470 mm (14.2 x 27.2 x 18.5 in)
- Weight: approx. 50 kg (110 lbs)

- Recommended air pressure: 4.5 bar (65 psi)
- Compressed air consumption at recommended air pressure: approx. 35 l/min (1.2 cfm)
- Air purity according to ISO 8573-1:2010
- 4-axes mechanism; positioning range (x/y/z): 80 x 60 x 40 mm (3.2 x 2.4 x 1.6 in)
- Built-in tank for cooling liquid in drawer shape with coarse filter mat, max. 4.2 l (4.4 qt)
- Complete housing of the working chamber; working chamber door with a safety contact and safety interlock
- Flow sensor for monitoring the cooling liquid circuit
- Working chamber lighting with different colours to indicate the machine state

Rotary axes

- Rotational range: +190 ° to -10 °
- Integrated blank holder for 3 blanks

Spindle

- Synchronous spindle rotating up to 60,000 RPM
- Nominal power under continuous load (S1): 300 W
- Nominal power under uninterrupted periodic operation (S6): 450 W
- Maximal Power (P_{max}): 600 W
- 4-fold bearing
- Hybrid ceramic ball bearing
- Pneumatic collet chuck ø 3 mm
- 8 cooling liquid nozzles
- Cone cleaning and sealing air

Tool magazines and tool change

- Automatic tool change
- 8 tool positions distributed on 2 exchangeable tool magazines
- Maximum tool length: 35 mm
- Automatic tool length measurement and tool breakage monitoring via measuring key
- Compressed air monitoring for the automatic tool change

Cooling lubricant

- Mix ratio with water: 1:19 (5% solution)
- pH value of 9.5

4 Installing the machine

4.1 Checking the scope of delivery

M Unpack the machine and ensure that you have received the following items:

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

- [1] 1 machine N4 Impression
- [2] 1 spindle service set
- [3] 1 power supply cord
- [4] 1 USB cable
- [5] 1 service unit for compressed air
- [6] 1 tank of pH testing strips
- [7] 1 bottle of cooling lubricant
- [8] 1 pneumatic hose
- [9] 1 torque wrench (1.8 Nm)
- [10] 2 keys for the emergency release of the working chamber door
- [11] 2 coarse filter mats (1 installed, 1 as a spare part)
- [12] 1 calibration set: 1 micrometer, 4 blanks, 2 pins (already screwed into a blank), 1 calibration tool
- [13] 1 measuring pin
- [14] 1 drill 2.8 mm for tool positions
- [15] 2 tool magazine inserts (as spare parts)

Without illustrations:


- This document
 - 1 supplement about removing the transport lock
 - 3 screws for the blank holder as spare parts
 - 1 measuring cup (250 ml)
- Please keep the machine packaging for possible further service shipments.

4.2 Choosing the installation site


- Choose the installation site according to the following criteria:
- Firm and even surface, has to carry the weight of the machine
 - Ambient temperature (storage / transport): between -20 °C (-4 °F) and 60 °C (140 °F)
 - Ambient temperature (operation): between 10 °C (50 °F) and 32 °C (90 °F)
 - Relative air moisture maximum 80 %, non-condensing
 - Machine location must be dust-free
 - Alternating current source with 100 – 240 V and 50/60 Hz
 - Incoming air pressure: maximum of 8 bar
 - Compressed air supply that meets the requirements of the machine (☞ chapter 3.7, page 10)
 - Sufficient space at the sides of the machine (☞ chapter 4.3, page 12)

4.3 Machine installation (scheme)

Connections

Distances to maintain

4.4 Installing the pneumatics

WARNING

Risk of injuries through leaking compressed air and lashing pneumatic hoses

Open or loose pneumatic connections can cause severe injuries.

- Make sure that **during installation** of the pneumatic hoses and of the service unit compressed air is not running through the hoses and connections.
- After installing the pneumatic hoses but **before** running compressed air through the hoses and connectors, check if the hoses are securely inserted into the correct connectors and are not damaged.
- Do not run compressed air through damaged hoses and connectors.

NOTICE

The spindle may suffer bearing damage and electrical damage if the compressed air is contaminated

The incoming compressed air must be **dry** and **oil-free** according to ISO 8573-1 because the service unit only serves as an **indicator** for contaminated air.

Air purity according to ISO 8573-1		
Solid particles	Class 3	Filtration degree better than 5 µm for solid particles
Water content	Class 4	Maximum pressure dew point +3 °C
Residual oil content	Class 3	Maximum oil content: 1 mg/m ³

- Make sure that the compressed air meets the mentioned requirements.
- Only connect the machine to the compressed air supply via the provided service unit.

The spindle requires compressed air for the following tasks:

- For the opening and closing of the collet chuck during tool change.
- For the sealing air which prevents foreign bodies from entering the spindle.

Air consumption of the machine:

- approx. 35 l/min at 4 bar (1.2 cfm at 60 psi)
- approx. 50 l/min at 8 bar (1.8 cfm at 120 psi)

4.4.1 Overview service unit

Via the service unit you connect the machine to your compressed air supply and regulate the incoming pressure for the machine.


FIG. 5: THE SERVICE UNIT: REGULATING AND CHECKING THE AIR PRESSURE

- [1] Rotary knob for pressure regulation
- [2] Pneumatic connection to the machine (ø 6 mm)
- [3] Manometer for monitoring the outgoing pressure
- [4] Water separator
- [5] Discharging screw
- [6] Standard pneumatic connection

4.4.2 Mounting the service unit on the machine

NOTICE

Failure of the water separator caused by a wrong alignment of the service unit

The service unit must **always** be mounted **in an upright position** because otherwise the water separator will not work.

- Mount the service unit in an upright position (↕ Fig. 5).

On the left side of the machine are two drillings which you can use to mount the service unit on the machine.

- M1. Remove the two blind screws in the drillings.
- M2. Mount the service units with the provided lens head screws in an upright position.


FIG. 6: MOUNTING THE SERVICE UNIT ON THE MACHINE

4.4.3 Installing the pneumatic hoses to the service unit


FIG. 7: INSTALLATION OF THE SERVICE UNIT (SCHEMA)

- M1. Shut off the external compressed air supply.
- M2. Via the provided pneumatic hose connect the right pneumatic connection of the service unit [2] with the pneumatic connection of the machine.
- M3. Connect the external compressed air supply with the standard pneumatic connection of the service unit [6].
- M4. Check the installation carefully for errors and damage. Do not conduct compressed air through damaged or loose hoses and connections!
- M5. Open the external compressed air supply.
- M6. Check the air pressure on the manometer [3]. If it does not lie between 4 bar and 8 bar, adjust it with the service unit (↪ chapter 4.4.4).

4.4.4 Setting the air pressure with the service unit

Setting the air pressure is only necessary if the air pressure shown by the manometer does not lie between 4 bar (60 psi) and 8 bar (120 psi).

- ★ The service unit is connected to your compressed air supply.

- M1. Pull the rotary knob on top of the service unit a little bit upwards.
- M2. Turn the rotary knob in the desired direction until the pressure lies between 4 bar and 8 bar – 4.5 bar (65 psi) are recommended:
 - Turning it towards „+“ you increase the pressure
 - Turning it towards „-“ you decrease it
- M3. Push the rotary knob down again.
 - The knob is locked and cannot be changed inadvertently.


FIG. 8: SETTING THE AIR PRESSURE

4.5 Filling the cooling liquid tank

- M1. If the working chamber door is open, close it.
- M2. Pull the cooling liquid tank out of the machine as far as it goes.

FIG. 9: PULLING OUT THE COOLING LIQUID TANK

- M3. Fill the cooling liquid tank with 4 litres of tap water.
- M4. Add 0.2 litres of cooling lubricant.
 - You have a 4.2 litre emulsion of water and cooling lubricant with a mixing ratio of approx. 5 %.



FIG. 10: MIXING WATER AND COOLING LUBRICANT

- M5. Take a pH testing strip from the provided tank and immerse it into the cooling liquid for a few seconds.

! Close the tank with the pH testing strips completely after usage. Otherwise the testing strips will become useless when exposed to air humidity for too long.

- M6. Compare the four coloured squares on the testing strip with the scale on the container. You can find the determined pH value in the column of the scale where the colours match the colours on the testing strip.

FIG. 11: COMPARISON OF pH SCALE WITH TESTING STRIP
FOR TECHNICAL REASONS COLOURS MAY DEVIATE

- [1] Determined pH value (exemplary value)
[2] Scale on container
[3] Testing strip

- M7. If the pH value lies beneath 9, add 50 millilitres cooling lubricant to the emulsion. Repeat steps M5 to M7, but no more than twice. After that dispose of the cooling liquid completely (↗ chapter 6.9, page 29) and start again with step M3.

- M8. If the pH value of the cooling liquid is correct, close the cooling liquid tank.

4.6 Establishing the electric connection

NOTICE

Damaging of the machine through heavy voltage fluctuations

Heavy voltage fluctuations can disrupt the control unit and can cause a failure of the system.

- Plug the machine's power cord in a dedicated circuit current or ensure that no devices are connected that can cause heavy voltage fluctuation when switched on.

NOTICE

Damaging of the machine if the transport lock and the CAM computer are installed

When you connect the machine to the electrical source and the CAM computer is connected, the machine starts referencing. During this process, the transport lock which is installed at delivery can damage the mechanics of the machine.

- Do not connect the machine to the electrical source if the CAM computer **and** the transport lock are installed.
- If the transport lock is installed, disconnect the USB connector between the machine and the CAM computer **before** connecting the machine to the electrical source.

- M1. Plug the provided power cord into the power connection at the connection panel of the machine.

- M2. Put the plug of the cord into a socket that is protected by a Residual Current Device / Ground Fault Circuit Interrupter.

4.7 Removing the transport lock

Before operating the machine for the first time, you must remove the transport lock. The transport lock prevents the spindle from getting damaged during transport.

- M1. Disconnect the USB connector between the machine and the CAM computer.

- M2. Connect the machine to the electrical source.

- M3. Switch on the machine.

- M4. Open the working chamber door.

- M5. Carefully pull the transport lock upwards and then out of the working chamber.

- M6. Clean the working chamber from parts of the transport lock that may have broken off.


FIG. 12: REMOVING THE TRANSPORT LOCK

4.8 Connecting the CAM computer

! The transport lock must not be installed when you connect the CAM computer to the machine (↗ chapter 4.7, page 15).

- M1. Switch on the machine.
- M2. Close the working chamber door.
- M3. Start the CAM computer.
- M4. Use the provided USB cable to connect a USB port of your computer or docking station with the connection panel of your machine.
- S5. Install the newest version of DentalCAM and DentalCNC that is released for the machine. For more information on this, read the documentation for the applications.
- S6.  Start DentalCNC and click on the depicted icon in the icon bar.
- S7.  Click on the depicted icon in the lower icon bar.
 - The application settings view displays.
- S8.  Click on the depicted icon beside the text Port number.
 - DentalCNC will try to establish a connection to the machine. If this is successful, the application will display the port number beside the icon and the machine references.
- !** The machine will not reference if the working chamber door is open.
- M9. If the working chamber door was open at step S8, close it. If the machine does not reference as a result, quit DentalCNC and restart the application.
- M10. If DentalCNC did not determine the port number, read how to determine the port number manually in the documentation for the application.

4.9 Testing the machine

After the first installation or a re-installation, for example after a transport, you should test the basic functions of the machine.

- ★ You installed the machine, the CAM computer and the manufacturing software completely. The machine is supplied with compressed air.
- M1. Start the machine.
 - The working chamber is illuminated in white.
- S2. Start DentalCNC.
- S3. Enter the USB port number in the application settings or have the program detect the port. For more information, read the documentation for DentalCNC.
 - The machine references. After that, the working chamber is illuminated in green. The machine is now in *default position*.
- S4. Switch to the machine control view.
- M5. Make sure that enough cooling liquid is in the cooling liquid tank.

NOTICE

Damaging of the pump when operated without cooling liquid

If the pump is operated without cooling liquid, the pump motor may overheat and get damaged.

- **Before** you activate the cooling liquid supply system, check if there is a sufficient amount of cooling liquid in the cooling liquid tank.
- If no cooling liquid is sprayed into the working chamber within a few seconds, deactivate the cooling liquid supply system **at once**.
- Do not use the machine for wet machining until the cooling liquid supply system is operational again.

- S6.  Activate the cooling liquid supply system by clicking on the depicted icon.
 - Cooling liquid sprays into the working chamber.
- S7. Use the depicted icons to verify that the compressed air and the cooling liquid flow rate are sufficient.
 - The icons display in blue when the compressed air and the cooling liquid flow rate are sufficient.

FIG. 13: LEFT: ICON FOR COMPRESSED AIR
RIGHT: ICON FOR FLOW RATE

- S8.  Move the spindle into the manual tool change position by clicking on the depicted icon.
- The spindle moves through the working chamber at constant speed.
- S9.  Move the spindle back to the default position by clicking on the depicted icon.
- The spindle moves to the default position at constant speed.
- S10.  Deactivate the cooling liquid supply system by clicking on the depicted icon.
- Cooling liquid no longer sprays into the working chamber.
- M11. Open and close the working chamber door.
- The working chamber door can be opened and closed easily.
- M12. If a result did not occur as described, check according to the error:
- The USB connection & USB driver installation ([↗](#) chapter 4.8 & documentation for the manufacturing software)
 - The compressed air supply ([↗](#) chapter 4.4)
 - The cooling liquid in the cooling liquid tank ([↗](#) chapter 4.5)
- M13. If you cannot solve a problem that occurred, contact customer service.

5 Running the machine

5.1 Switching on the machine

- ★ The machine is installed completely and is supplied with compressed air. The working chamber door is closed.
- M1. Switch on the machine and the CAM computer.
- The working chamber is illuminated in white.
- S2. Start DentalCNC.
- The machine references. After that, the working chamber is illuminated in green. The machine is now ready for use.

5.2 Overview job execution

NOTICE

Damaging of the machine when using damaged tools or blanks

If tools or blanks are damaged, parts can break off and damage the machine during job execution.

- Check the blanks and tools **thoroughly** for damage before the job execution.

Carrying out a typical job with the N4 Impression is as follows:

1. Select the blanks and the jobs in DentalCAM
2. Insert the tools into the machine
3. Mount the blanks into the machine
4. Execute the job in DentalCNC

The steps in the manufacturing software are described in the corresponding documentation. The following text describes how to perform manual work with the machine.

5.3 Inserting / changing tools

NOTICE

Damaging of the spindle or the tool positions if you use improper tools

Improper tools can damage the collet chuck of the spindle and /or the tool positions.

- Only use tools with a sufficiently large chamfer at the tool shank.
- Install a retaining ring as a stop ring according to DIN 471-A3.
- Only insert tools with a maximum diameter of 3 mm at the thickest part into the collet chuck.
- Only insert tools with a maximum cutting edge diameter of 2.6 mm into the tool magazines.

 We recommend original tools as they are designed especially for the designated jobs.

You can insert up to 8 tools into the tool magazines in the working chamber of the machine. The machine can change tools automatically during machining so that it executes jobs without your interference.

The 8 tool positions are located on two extractable tool magazines inside the N4 Impression.


FIG. 14: TOOL MAGAZINES IN THE WORKING CHAMBER OF THE MACHINE
TOP MARKING: TOOL MAGAZINE A
BOTTOM MARKING: TOOL MAGAZINE B

You can equip the tool magazines in two ways:

- Via the spindle – you insert a tool into the collet chuck and the spindle deposits the tool in the tool magazine afterwards. This function is described in the documentation of the manufacturing software.
- By inserting the tools into the tool magazines manually. This is described below.

You insert the tools manually as follows:

- ★ The machine is in default position.
No items are mounted in the blank holder.

No tool is inserted in the collet chuck.

- M1. Open the working chamber door of the machine and pull out the tool magazines to the right.



FIG. 15: PULLING OUT THE TOOL MAGAZINES

- S2. Start DentalCNC. Call up the tool administration view.
- S3. Select the tool magazine with the tools that you want to insert into the machine or create a suitable magazine (e.g. for the next job).
- M4. Insert the tools into the tool positions on the *upper side* of the tool magazines:
- The positions of the tools in the blocks must match the tool positions in DentalCNC ( Fig. 16).
 - Insert the tools straight into the tool positions and push them in until the ring touches the rubber ( Fig. 17).
-  If the positions of the tools in the tool magazines do not match the tool positions in DentalCNC, the machine will use the wrong tools during job execution and the job result will become unusable.


FIG. 16: COMPARISON OF THE CNC SOFTWARE AND TOOL MAGAZINES

- [1] Tool positions in DentalCNC
- [2] Tool positions on the block
- [3] Tool magazine insert


FIG. 17: INSERTING TOOL STRAIGHT INTO THE TOOL POSITION; RING AT THE TOOL ENLARGED

NOTICE

Damaging of the tool magazine holder or the tool magazines through wrong insertion or insertion with force

Each tool magazine only fits into its designated place in the tool magazine holder. The blocks are mechanically protected so that the user cannot completely insert them into the wrong place.

- Only insert a block into the place from which you removed it.
- Do not push the tool magazines into the tool magazine holder with force.

M5. Re-insert the tool magazines into the tool magazine holder and ensure that they are properly seated.

- You can execute the jobs for the inserted tools.

5.4 Mounting and removing blocks

You can mount up to 3 blocks into the machine and have them processed in one job.

⚠ The blocks for one job need to be of the same type (material and dimensions).

- M1. Open the working chamber door of the machine.
- M2. Loosen the screws of the blank holder with a hexagonal socket screwdriver and remove already mounted items (if any). You do not need to remove the screws from the blank holder completely.
- M3. Insert the block mandrels into the openings in the blank holder and push the blocks into the blank holder until they are properly seated.

- M4. Immobilize the blocks with the screws which you loosened in step M2. Use the provided torque wrench (1.8 Nm) for this and tighten the screws until the torque wrench triggers.


FIG. 18: INSERTING THE BLOCKS

5.5 Using an optional abutment holder

With the optional abutment holders your machine can process abutments with prefabricated connection geometries.

- i** Prefabricated abutments usually require manufacturer specific abutment holders. Different abutment holder models can be obtained from customer service.

5.5.1 Installing the PreFace® abutment holder

- S1.  Move the spindle back to the default position by clicking on the depicted icon.
- M2. Open the working chamber door of the machine and remove all items from the blank holder.
- M3. Switch off the machine.
- You can manually rotate the blank holder.
- M4. Rotate the blank holder into a vertical position (↗ Fig. 19).
- M5. Place the abutment holder onto the blank holder and screw the long screws of the abutment holder into the blank holder manually.


FIG. 19: SCREWING IN THE ABUTMENT HOLDER

- M6. Turn the blank holder clockwise up to the stop.
- M7. Screw down the 2 screws with the provided ratchet.


FIG. 20: IMMOBILIZING THE ABUTMENT HOLDER

- M8. Close the working chamber door of the machine.
- M9. Switch on the machine.

5.5.2 Demounting the PreFace® abutment holder

- M1. Remove inserted prefabricated abutments from the abutment holder (if any).
- M2. Loosen the screws of the abutment holder (↖ Fig. 20) and remove it.

5.6 Mounting and removing PreFace® abutments

- M1. Loosen the screw at the PreFace® abutment holder with a hexagonal socket screwdriver. You do not need to remove the screw completely.
- M2. Carefully push the loosened screw of the PreFace® abutment holder downwards into the holder.


FIG. 21: OPENING THE PREFACE® ABUTMENT HOLDER

M3. Remove previously mounted PreFace® abutments.

NOTICE

Bad machining results caused by incorrect scanning of prefabricated abutments

If moisture or dirt covers the prefabricated abutments, the machine will scan them incorrectly and the machining result will become unusable.

- Wipe over the bellow before you insert prefabricated abutments so that no cooling liquid drips onto them.
- Clean the prefabricated abutments with a dry, lint-free cloth before processing.

M4. Push the PreFace® abutments into the PreFace® abutment holder until they are properly seated.

M5. Retighten the screw at the PreFace® abutment holder.

- The PreFace® abutments are immobilized.


FIG. 22: INSERTING THE PREFACE® ABUTMENTS

5.7 Executing jobs

M1. Make sure that the following requirements are met before job execution:

- You created the job that you want to execute in DentalCAM and transferred it to DentalCNC.
- All required tools are assigned to the job in DentalCNC.
- The required tools are in the right positions in the tool magazines of the machine.
- All required blanks of the correct material are inserted.
- The cooling liquid tank is filled with a sufficient amount of cooling liquid.

! Do not move the machine during job execution because otherwise the results may become imprecise.

M2. Close the working chamber door and the cooling liquid tank.

S3. Start the job execution in DentalCNC.

- The machine begins with the job execution. The working chamber is illuminated in blue.

5.8 Interrupting the job execution

S  Click on the depicted icon to stop the job execution. The icon only displays during job operation.

5.9 Operation interruption or abortion of job execution

An operation interruption can occur because the compressed air supply or the vacuum in the working chamber is not sufficient. An *interrupted* job will normally be resumed automatically after the error is corrected.

The job execution is *aborted* in the following cases:

- In case of a machine malfunction
- In case of a tool breakage
- In case of a power failure

If a job was *aborted*, you have to execute it again in DentalCNC.

5.9.1 How to proceed in case of an operation interruption

If the job execution was interrupted, DentalCNC displays a corresponding message.

M1. If DentalCNC prompts that the compressed air is insufficient, check the following:

- The manometer of the service unit
- The installation of the pneumatic hoses
- Your compressor

M2. If DentalCNC prompts that the flow rate of the cooling liquid is too low, ensure that the cooling liquid tank contains a sufficient amount of cooling liquid.

 You will find additional information in the chapter on troubleshooting ([↗ chapter 8, page 36](#))

5.9.2 How to proceed in case of a machine malfunction

A machine malfunction is recognized by the internal control unit in case of a critical event. The working chamber will be illuminated in red. DentalCNC displays the error message and error code that was sent by the control unit.

M1. Disconnect the machine from the electrical source and prevent it from being restarted.

M2. Record the error message and error code that is displayed in DentalCNC.

M3. Contact customer service. Have the error message and error code from DentalCNC readily available.

M4. To remove a blank from the working chamber, unlock the working chamber door with the key for emergency release ([↗ chapter 2.2.1, page 6](#)).

5.9.3 How to proceed in case of a tool breakage

If a tool breaks during machining, the machine will not recognize this immediately and the spindle will continue until the next tool change. Only upon measuring the broken tool will the job execution be cancelled and DentalCNC will display a corresponding message.

A tool breakage can be caused by the following:

- The tool was damaged or worn.
- The tool was put into the wrong tool position or was inserted into the spindle manually at the wrong time. As a consequence, it was not suitable for the processing step.
- The distribution of the objects in the blank (“nesting”) was not sufficiently adapted to the material.

M1. Open the working chamber door.

M2. Remove all parts of the broken tool from the working chamber.

M3. If the spindle picked up the tool from the tool magazines, check if the tool was inserted into the correct position. Insert a spare tool into the correct position in the tool magazines ([↗ chapter 5.3, page 18](#)).

M4. If you inserted the tool into the collet chuck manually, check if the broken tool corresponds to the initial tool type in which DentalCNC prompted. Have a correct spare tool ready.

M5. Close the working chamber door and execute the job again.

 If tools break regularly, you can find additional information in the chapter on troubleshooting ([↗ chapter 8, page 36](#))

5.9.4 Procedure in case of a power failure

As long as the machine is not powered, the working chamber door will remain locked.

- After a short power failure turn on the machine again and follow the on-screen prompts so that the machine can reference.
- In case of a longer power failure or a machine malfunction use the provided key for the emergency release of the working chamber door.

5.10 Emergency release of the working chamber door

- Only use the emergency release of the working chamber door to remove a blank from the working chamber if the working chamber door is permanently locked.

WARNING

Crushing hazard and cutting injuries while running the machine with the working chamber door open

If the safety lock is not reactivated after the emergency release, the user will be able to operate the machine with the working chamber door open. In this case the working chamber door will not protect the user from bruises caused by machine movements and from cutting injuries caused by the rotating tool.

- Activate the safety interlock **immediately** after completing your work in the working chamber.
 - **Never** operate the machine with the working chamber door open.
 - Perform the emergency release only if you are authorised to do so and if you have received training on its usage.
 - Store the key for the emergency release in a place where only authorised persons can access it.
 - Report **every** emergency release to customer service.
-

CAUTION

Cutting injuries when touching a rotating tool

If a power failure or a machine malfunction occurs during processing, the spindle including the inserted tool keeps rotating. If you touch the rotating tool, you will suffer from cutting injuries.

- Wait until the spindle has stopped rotating before carrying out an emergency release.

There is an opening for the emergency release at the left side of the machine.


FIG. 23: OPENING FOR THE EMERGENCY RELEASE

- M1. Switch off the machine at the main power switch and disconnect the machine from any electrical source and the compressed air supply.
 - M2. Remove the protective cap from the emergency release opening at the left side of the machine (☞ Fig. 24, top).
 - M3. Insert the key for the emergency release of the working chamber door through the opening at the left side of the machine. (☞ Fig. 24, middle).
 - M4. Deactivate the safety interlock of the working chamber door by turning the key counter-clockwise by 90° (☞ Fig. 24, bottom).
- You can open the working chamber door.
- M5. Carry out your work in the working chamber.
 - M6. Close the working chamber door.
 - M7. Reactivate the safety interlock of the working chamber door by turning the key for the emergency release clockwise by 90° (☞ Fig. 24, bottom).

- M8. Check if the working chamber door actually is locked.
 - If you can still open the working chamber door, repeat the steps M6–M8.


FIG. 24: USING THE KEY FOR THE EMERGENCY RELEASE

TOP: REMOVE CAP
MIDDLE: INSERT KEY
BOTTOM: TURN THE KEY

- M9. Close the opening for the emergency release again with the protective cap.

6 Machine maintenance and cleaning

6.1 Definition of wear parts

The machine and the extra equipment are warranted for a period of 24 months or 2000 operating hours, whatever comes first. The warranty covers defects in materials or fabrication as long as the regulations for using the machine in all documents are followed.

Of course, the warranty also covers wear parts as long as their failure cannot be attributed to the function-related abrasion. The wear parts that are mentioned below can already wear down within the warranty period due to their normal function. The average useful life of the wear parts can be seen in the following table.


Use these values to determine operating costs, to plan your spare part stock as well as to create individual maintenance and service plans.

6.2 Maintenance table

Task	Recommended interval	Procedure / Utensils	Illustration spare part
Cleaning the working chamber & the measuring key (☞ chapter 6.4, page 25)	1 x per day, more often if the chamber is very dirty	Wet vacuum cleaner, brush, wet cloth, never use compressed air	
Cleaning the collet chuck (☞ chapter 6.5, page 26)	1 x per week, in case of rotational imperfections, in case of bad machining results	Spindle service set	
Checking the hoses, cables & connections (☞ chapter 6.6, page 27)	1 x per week	Visual inspection	
Changing the cooling liquid, cleaning the cooling liquid tank and the coarse filter mat (☞ chapter 6.9, page 29)	10-20 operating hours*	Cloth, water & cooling lubricant, pH testing strips, possibly a bucket	
Cleaning the housing (☞ chapter 6.11, page 30)	As necessary	Cloth, water, possibly a mild cleaner	
Service unit: a) Checking for contamination b) Cleaning / exchanging the filter cartridge (☞ chapter 6.12, page 31)	a) Necessarily every day b) With visible pollution or every 2 years*	Possibly a new filter cartridge	
Exchanging the main fuse (☞ chapter 6.13, page 31)	In case of a defect fuse in the power supply of the machine	Replacement fuse of the type T6,3A L250V	
Calibrating the axes (☞ chapter 6.14, page 32)	Only if necessary	Calibration specimen, micrometre	
Inspection by service technician (☞ chapter 6.3, page 25)	After 2000 operating hours* or every 2 years*		

* These are recommendation guidelines. Depending on the processing material and how well the machine is cleaned, these values may differ.

Wear parts

Exchanging wear parts	Average useful life	Procedure / Utensils	Illustration spare part
Replacing the coarse filter mat (☞ chapter 6.9, page 29)	100 operating hours*	New Coarse filter mat	
Exchanging the tool magazine inserts (☞ chapter 6.15, page 34)	500 operating hours*	Replacement inserts, drill for drilling the tool positions, screwdriver	
Exchanging the collet chuck	1000 operating hours*		
Exchanging the pump (☞ chapter 6.10, page 30)	250 operating hours*		
Exchanging the spindle bearings (requires spindle exchange by service technician)	2000 operating hours* or every 2 years*	Performed by service technician	

* These are recommendations as a guidance. Depending on the processing material and pollution of the machine these values can differ.

6.3 Inspection

We recommend having a service technician perform an inspection regularly.

- In addition, have an inspection carried out whenever the machine is being repaired.

6.4 Cleaning the working chamber and measuring key

NOTICE

Damaging of the linear guides or the spindle when cleaning with compressed air

If you clean the working chamber with compressed air, material chips can reach the linear guides or the spindle bearings.

- Clean the interior of the machine **exclusively** with a wet vacuum cleaner, a brush and a cloth and **never** with compressed air.

M1. Have ready:

- A wet vacuum cleaner
- A wet cloth
- A *mild* cleaner (optional)
- A wet brush

M2. Turn on the machine.

S3. Start DentalCNC.

S4.  Activate the sealing air by clicking on the depicted icon.

- The sealing air is audibly blown into the working chamber.

M5. Open the working chamber door.

M6. Suck up the coarse dirt in the working chamber with the wet vacuum cleaner.

M7. Clean the whole interior thoroughly with the wet cloth. Use a *mild* cleaner if necessary.

M8. Thoroughly clean the blank holder including all screws, pins and clamping mechanisms as well as all openings from the processing dust.

M9. Clean the measuring key between the tool magazines with a brush.


FIG. 25: MEASURING KEY BETWEEN THE TOOL MAGAZINES (MARKED ORANGE)

M10. Close the working chamber door.

S11.  Deactivate the sealing air by clicking on the depicted icon.

- The sound caused by the sealing air subsides.

6.5 Cleaning the collet chuck

NOTICE

Damaging of the spindle when cleaning with compressed air

If you clean the collet chuck with compressed air, the spindle bearings may get damaged.

- Clean the collet chuck **exclusively** with the appropriate service set.


FIG. 26: COMPONENTS OF THE SPINDLE SERVICE SET

- [1] Collet grease
- [2] Cleaning brush
- [3] Knurled nut
- [4] Cleaning cone

- M1. Close the working chamber door of the machine.
- S2. Open the collet chuck with the “Release tool” function of DentalCNC and remove the tool from the spindle if necessary.
- M3. Open the working chamber door and remove the tool from the spindle if necessary.
- M4. Insert the measuring pin into the collet chuck and keep holding it with one hand. With the other hand put the knurled nut to the spindle.


FIG. 27: UNSCREWING THE COLLET CHUCK

- M5. Loosen the collet chuck with the knurled nut and unscrew the collet chuck with your hand or with the knurled nut.
- M6. Remove the measuring pin from the collet chuck. Put it aside within reach together with the knurled nut.
- M7. Clean the inner cone of the spindle with the cleaning cone of the service set.


FIG. 28: CLEANING THE INNER CONE OF THE SPINDLE

- M8. Clean the collet chuck with the brush of the service set. Make sure that the 3 buffers in the collet chuck do not fall out during cleaning.


FIG. 29: CLEANING THE COLLET CHUCK

- [1] Longitudinal slot
[2] Buffer
[3] Cleaning brush

NOTICE

Damaging of the spindle when using the wrong or wrongly applied grease

If you use grease that is unsuitable for the spindle or if grease gets into the longitudinal slots of the collet chuck, the spindle can get damaged.

- Make sure that no grease gets into the longitudinal slots of the collet chuck.
- Only use a very small, about pinhead-size amount of the grease.
- Only use the provided collet grease of the service set.

- M9. Put a small amount of the collet grease on the index finger and smear it with the thumb.
- M10. Apply the smeared collet grease to the flanks of the collet chuck.


FIG. 30: GREASING THE COLLET CHUCK; SURFACE TO GREASE MARKED BLUE; SLOT THAT MAY NOT BE GREASED MARKED RED

- M11. Insert the measuring pin into the collet chuck with one hand and keep holding it. Screw the collet chuck with the knurled nut in your other hand *tightly* into the spindle.

! In any case, turn the knurled nut up to the stop so that the collet chuck is properly seated in the spindle. Otherwise rotational imperfections can occur during operation which will worsen your processing results.

- M12. Remove the measuring pin from the collet chuck and store it together with the other components of the spindle service set.

- M13. Close the working chamber door.

6.6 Checking the hoses, cables and connections

- M1. Disconnect the machine from the electrical source and the compressed air supply.
- M2. Check the hoses and cables that are installed to the machine.
- M3. In case of damages prevent the machine from being restarted. Replace the damaged cables and hoses immediately with original replacement cables. Contact customer service.
- M4. Check all connections at the connection panel for loose connections.
- M5. Re-insert loose cables and hoses tightly into the corresponding connections. In case of damaged connections prevent the machine from being restarted immediately. Contact customer service.

6.7 Removing the cooling liquid tank completely from the machine

You can remove the cooling liquid tank completely from the machine, e.g. to clean it thoroughly on a bench or to pour the cooling liquid into a container. When exchanging the pump, you also have to remove the cooling liquid tank.

M1. Close the working chamber door if necessary.

M2. Open the cooling liquid tank. Pull it out up to the stop.

M3. Pull the black locking lever at the left side of the machine upwards and keep it in this position.

M4. Push the black locking lever at the right side of the machine downwards and keep it in this position.

M5. Pull the cooling liquid tank carefully out of the machine and let go of the locking levers.


FIG. 31: LOCKING LEVERS AT THE SIDES OF THE COOLING LIQUID TANK

[1] Locking lever at the left side

[2] Locking lever at the right side

6.8 Re-inserting the cooling liquid tank


FIG. 32: INSERTING THE COOLING LIQUID TANK INTO THE MACHINE

M1. Close the working chamber door.

M2. Insert the guide rails at the left and the right side of the cooling liquid tank into the guide rails of the machine. Afterwards push the tank carefully into the machine.

6.9 Exchanging the cooling liquid / cleaning the cooling liquid tank

NOTICE

Increased wear of the pump if the cooling liquid tank is not sufficiently cleaned

If the cooling liquid tank is not sufficiently cleaned before filling in the new cooling liquid, the useful life of the pump can decrease significantly.

➤ **Every time** that you exchange the cooling liquid, carry out **all** steps in this chapter in the given order.

- M1. *Before* exchanging the cooling liquid and /or cleaning the cooling liquid tank, clean the working chamber and measuring key.
- M2. Have ready:
 - A container for collecting the cooling liquid
 - A cloth for cleaning
 - Tap water
 - A wet vacuum cleaner if necessary
- M3. Close the working chamber door.
- M4. Switch off the machine.
- M5. Carefully remove the cooling liquid tank completely from the machine (☞ chapter 6.7, page 28).
- M6. Completely empty the cooling liquid tank, filling its content into the collecting container.
- ! For disposing the cooling liquid read the corresponding chapter (☞ chapter 7.1, page 35).
- M7. Pull the coarse filter mat upwards and out of its holder.


FIG. 33: PULLING THE COARSE FILTER MAT OUT OF ITS HOLDER

- M8. Clean the cooling liquid tank thoroughly with a cloth and a lot of tap water.
- M9. If you want to reuse the coarse filter mat, clean it thoroughly with a lot of tap water.

- M10. Remove the residual dirt in the cooling liquid tank with a wet vacuum cleaner if necessary.
- M11. Insert the cleaned mat or a new mat into the holder.
- M12. Fill tap water without cooling lubricant into the cooling liquid tank until the pump is completely submerged.
- M13. Re-insert the cooling liquid tank into the machine (☞ chapter 6.8, page 28).
- M14. Switch on the machine and start DentalCNC.
- S15.  Activate the cooling liquid supply system by clicking on the depicted icon.
 - Cooling liquid is spraying in the working chamber.
- M16. Check the purity of the water through the view window to the working chamber. Do not continue as long as the water is not clear.
- S17.  Deactivate the cooling liquid supply system by clicking on the depicted icon.
 - Cooling liquid is no longer spraying into the working chamber.
- M18. Switch off the machine.
- M19. Dispose of the tap water in the cooling liquid tank and re-insert the cooling liquid tank into the machine.
- M20. Fill water and cooling lubricant in the cooling liquid tank as described in the corresponding chapter (☞ chapter 4.5, page 14).
- M21. Close the cooling liquid tank.

6.10 Exchanging the pump


FIG. 34: SPARE PART PUMP

- [1] Angle plate
- [2] Angle piece, permanently attached to pump
- [3] Pump
- [4] Contact plug

- M1. Remove the cooling liquid tank completely (↗ chapter 6.7, ↗ page 28).
- M2. Empty and clean the cooling liquid tank completely (↗ chapter 6.9).
- M3. Put the cooling liquid tank onto a firm and even surface with sufficient lighting.
- M4. Pull the pump filter from the defective pump.


FIG. 35: REMOVE PUMP FILTER

- M5. Clean the pump filter and push it onto the new pump up to the stop.
- M6. Loosen the 4 screws of the angle plate and the contact plug at the backside of the cooling liquid tank.


FIG. 36: BACKSIDE OF THE COOLING LIQUID TANK

- M7. Pull the cooling liquid hose out of the angle piece at the pump and let water that may leak flow into the cooling liquid tank.
- M8. Remove the pump including the angle plate and the contact plug.
- M9. Insert the cooling liquid hose into the angle piece of the new pump until it is properly seated.
- M10. Screw the angle plate of the new pump and the contact plug with the 4 screws to the backside of the cooling liquid tank.

6.11 Cleaning the housing

NOTICE

Damaging of the housing when choosing a non-recommended cleaner

If you use an unsuitable cleaner for cleaning the machine housing, the surface or the adhesive foil can get damaged.

- Avoid contact of the powder coating with alkaline or acid substances. Especially metallic powders show a high sensitive reaction.
- Take care of the glued symbols so that they do not peel off. The adhesive foil is especially sensitive to rubbing and strong detergents.
- If the usage of a special cleaner is necessary to eliminate certain dirt, we recommend checking the suitability of the cleaner at a hidden place of the part first.

- M1. First, clean the powder coated surface with a soft, dry cloth.
- M2. If some dirt cannot be removed this way, moisten the cloth. Use a pH neutral cleaner if necessary.

6.12 Checking the service unit

NOTICE

Damaging of the machine when compressed air is contaminated

Compressed air that does not fulfil the guidelines for purity according to ISO 8573-1 can damage the machine (↗ page 13).

- Check the water separator of the service unit daily for contamination.
- Never use the machine if there is water, oil or solid particles in the water separator.

6.12.1 Checking the water separator for condensate

Condensate in the separator usually points to compressed air not being dry enough.

- M1. Check if water, oil or solid particles piled up in the water separator.
 - In this case switch the machine off immediately and proceed as follows:
- M2. Check the compressed air supply and make sure that the compressed air fulfils the requirements for air purity according to ISO 8573-1 (↗ chapter 4.4.1, page 13). Do not use the machine until the compressed air fulfils this requirement!
- M3. Drain the water separator by turning the discharging screw counter-clockwise.
 - The condensate is blown out downwards under pressure.
- M4. Close the discharging screw again by turning it clockwise.


FIG. 37: OPENING/CLOSING THE DISCHARGING SCREW AT THE SERVICE UNIT

6.12.2 Exchanging / cleaning the contaminated filter cartridge

You have to clean or exchange the filter cartridge in the water separator in case of strong contamination.

- ❗ A strongly contaminated cartridge can lead to a pressure loss.

If the compressed air fulfils the requirements for air purity according to ISO 8573-1, the filter cartridge usually does not have to be changed.

- If the filter cartridge is contaminated, check the purity of your compressed air (↗ page 13).

You exchange or clean the filter cartridge as follows:

- M1. Disconnect the machine from the compressed air supply.
- M2. Unscrew the bowl of the water separator.
- M3. Untighten the filter screw below the filter cartridge.
- M4. Pull out the filter cartridge and clean it if necessary.
 - ❗ A new filter cartridge is available as spare part from customer service.
- M5. Insert the new or cleaned filter cartridge and reassemble the water separator.

6.13 Exchanging the main fuse

- ❗ Use only a fuse of the type T6,3A L250V as replacement fuse.

The internal power supply of the N4 Impression has a main fuse that is accessible from outside and can be replaced if necessary.

- ❗ A new main fuse is available as spare part from customer service.

- M1. Turn off the machine at the main power switch and disconnect the machine from any electrical source and the compressed air supply.
- M2. Remove the power cord from the connection panel.
- M3. Remove the cover of the fuse.


FIG. 38: COVER OF THE FUSE (MARKED ORANGE)

- M4. Remove the defective fuse and replace it with a new fuse of the type T6,3A L250V.

If you do not have a replacement fuse ready, take the replacement fuse from the right side of the fuse cover and put it into the left side of the fuse cover.

- M5. Remount the fuse cover.

6.14 Calibrating the axes

NOTICE

Deterioration of machining results caused by an incorrect calibration

The machine is already calibrated at delivery. As long as your machining results are without flaw, a new calibration is not necessary. A calibration takes much time and can deteriorate the machining results in the worst case if it is not done correctly.

- In case of inaccurate machining results, try adjusting the working conditions first: check the fixation of the blank, the state of the tool or of the processing material.
- Before calibrating the machine, contact customer service.
- Be very careful when measuring and entering data during calibration. When in doubt, stop the calibration.

i The documentation for the manufacturing software contains all information on calibrating the machine. Therefore, you only find information specific to the N4 Impression in this document.

Mounting a calibration / test specimen in the machine

Your N4 Impression is delivered with a calibration set. It contains the following parts:

- Calibration blanks which you use to mill calibration and test specimens
- Pins for immobilizing the blanks in the blank holder
- A tool for milling the calibration or test specimen
- A micrometre for measuring the machining precision


FIG. 39: BLANK AND PINS FROM THE CALIBRATION SET

- [1] Blank
- [2] Recess for pin
- [3] Pin for chucking the blank into the blank holder
- [4] Opening for the screw used to fixate the pin

i At delivery the pins are already screwed into one of the blanks. Use this blank first and, in this case, skip the first two steps.

NOTICE

Tool breakage or incorrect results caused by loose screws

If the blank is not properly seated in the blank holder, the measurement results may become incorrect or the tool may break when calibrating.

- Screw the pins tightly into the blanks.
 - Screw the blanks tightly into the blank holder.
- ★ The machine is in default position. There are no blanks mounted onto the blank holder.
- M1. Insert the pins into the recesses [2] in the blank. Make sure that the openings for the screws [4] are visible after the insertion.
- M2. Insert the screws into the openings of the pins and screw the pins tightly into the blank.


FIG. 40: INSERTING PINS AND SCREWING THEM IN TIGHTLY

- M3. Using the pins, mount the blank onto the blank holder in the working chamber of the machine.


FIG. 41: MOUNTING A BLANK FROM THE CALIBRATION SET INTO THE BLANK HOLDER

- S4. Mill the calibration/test specimen as described in the documentation for the manufacturing software.

! Calibration/test specimen are always machined dry. This does not indicate that the machine does not work correctly.

- M5. Remove the calibration/test specimen and the calibration tool from the machine.
 M6. Untighten the screws and remove the pins from the blank. Store the pins, screws and the calibration tool together for further usage.
 M7. Break the blank apart to measure it (see below) and follow the steps in the documentation for the CNC software.

Breaking the calibration/test specimen apart to measure it

After milling the calibration/test specimen looks as follows:


FIG. 42: A MILLED CALIBRATION/TEST SPECIMEN FOR THE N4 IMPRESSION

- M Break the blank into the following three pieces by pushing out the inner figure respectively the cube in the centre with your thumb:


FIG. 43: THE INDIVIDUAL PARTS OF THE TEST AND CALIBRATION SPECIMEN

- [1] Cube in the centre
- [2] Inner figure
- [3] Outer frame

6.15 Replacing the tool magazine inserts

When the tool magazine inserts are worn, they should be replaced. Afterwards you must drill the tool positions into the new inserts with the machine.

- Your machine comes with tool magazine inserts as spare parts and with the drill tool.
 - Additional inserts and drill tools are available via customer service.
 - Tool magazines are available from customer service as well.
- Always replace the inserts for *both* blocks together even if only one of the inserts is worn. The machine will always drill the tool positions for *both* inserts!
- i** The documentation for the manufacturing software contains step-by-step instructions for drilling the inserts. Below you find a description of how to replace the tool magazine inserts in the machine.

- M1. Have 2 tool magazine inserts ready as spare parts.
- M2. Open the working chamber door.
- M3. Pull one tool magazine out of the tool magazine holder and remove all tools.
- M4. Unscrew the 4 screws on the upper side of the block and lift the cover.


FIG. 44: UNTIGHTENING THE SCREWS


FIG. 45: LIFTING THE COVER

- M5. Remove the existing, worn tool magazine insert and replace it with a new one.


FIG. 46: NEW TOOL MAGAZINE INSERT (ON THE LEFT, WITHOUT DRILLINGS) REPLACING THE OLD INSERT (ON THE RIGHT, WITH DRILLINGS)

- M6. Put the cover back onto the tool magazine and fasten it with the screws that you unscrewed in step M3 (↻ Fig. 47).
- M7. Insert the tool magazine into the tool magazine holder.
- M8. Repeat the steps M3 to M7 with the other tool magazine.
- M9. Close the working chamber door.
- S10. Drill the tool positions into the new inserts.


FIG. 47: ASSEMBLING A TOOL MAGAZINE

7 Disposal

7.1 Disposing of the cooling liquid

- The product contains neither PCB nor other chloroparaffins.
 - Waste code number: 120109, machining emulsions and solution free of halogens
 - Burden of proof: yes
- Avoid entry of new or used machining emulsion into the soil, in waters or into the canalisation.
 - Recycle or incinerate the machining emulsion at authorised facilities. For this, transfer it to approved disposal companies.
 - Keep a reference sample of the disposal product for at least 6 months.
 - Abide the national and local laws of the disposal location in any case.

7.2 Disposing of the machine

We will dispose of the machine at no cost. The owner will bear the costs for disassembly, packaging and transport.

- Before sending in the machine for disposal, contact customer service
- If you dispose the machine yourself, abide the national and local laws of the disposal location.

8 Troubleshooting

DANGER

Danger to life due to improper troubleshooting

In case of improper troubleshooting, you may suffer from severe or deadly injuries caused by, for example, electric shocks, axis and spindle movements, loose pneumatic components and pneumatic hoses. In addition, your machine may get damaged.

- Do not remove the housing of the machine.
- **Before** you check or run pneumatic hoses, close the external compressed air supply valve.
- **Before** you check or run electric cables or power connectors, disconnect the machine from the electrical source and prevent it from being restarted.
- Do not perform any troubleshooting while the machine is operating.
- Wear safety gloves throughout the troubleshooting process.
- If you are unsure of how to perform certain steps during troubleshooting or cannot solve the problems, abort the troubleshooting and contact customer service.

Additional symbols in this chapter

- Problem / question
- Solution / answer

Frequently Asked Questions (FAQ)

- I cannot open the working chamber door.**
- Is the machine operating? While the axes are moving, the working chamber door stays locked.
- M Wait until the machine has finished and the working chamber lighting shines green.
- Without electricity, the working chamber door remains locked.**
- Has a power failure occurred at the installation site of the machine?
- M Depending on the duration of the power failure, restart the machine or use the key for the emergency release.
- chapter 5.9.4, page 22
- Is electricity available at the installation site of the machine?
- M1. If you have not removed the transport lock from the working chamber yet, disconnect the machine from the CAM computer.
- M2. Connect the machine to the electrical source and switch it on via the main power switch.

- M3. If the working chamber lighting does not illuminate, check if the power cable is properly seated in the machine and is connected to the power source.
- M4. Try connecting the machine to a different socket.
- chapter 4.6, page 15

I have installed all components, started DentalCNC but the machine does not reference.

- Is the working chamber door open? The machine does not reference with the working chamber door open.
- M Close the working chamber door.
- Is the USB cable connected properly?
- M Check if the USB cable is properly seated in the connector and is undamaged. If possible, use the provided cable.
- Is the correct port specified in DentalCNC? Without the correct port, the application will not establish a connection with the machine.
- S1. Start DentalCNC.
- S2. Try to have the application determine the correct port in the application settings.
- chapter 4.8, page 16
- S3. If the application cannot determine the correct port, enter the port manually.
- Documentation for the manufacturing software

My machine does not execute any jobs although DentalCNC displays that there is a connection between the CAM computer and the machine.

- Is the working chamber door open? The machine does not move the axes if the working chamber door is open.
- M Close the working chamber door.
- chapter 3.2, page 8

The machining results are not satisfactory and / or tools keep breaking.

- Do the tool positions in DentalCNC correspond to the tools in the tool magazines of the machine? If not, the machine uses the wrong tools during job execution.
- S1. Call up the job execution view in DentalCNC.
- S2. Compare the tool positions in the application to the tools in the tool magazines.
- M3. Replace the wrong tools in the tool magazines with the correct ones.
- chapter 5.3, page 18

 Is the blank attached properly?

M Ensure that all screws, pins and clamping mechanisms fixate the blank firmly.

 chapter 5.4, page 19

 Are the tool magazine inserts worn?

M Exchange the tool magazine inserts against new ones.

 chapter 6.15, page 34

 Are screws, pins and clamping mechanisms as well as the corresponding openings contaminated by processing dust?

M Clean the mentioned components thoroughly.

 chapter 6.4, page 25

 Are the tools worn?

M Check the tools in the tool magazines and exchange the worn tools with new ones. Also check the tool life values which are shown in DentalCNC.

 Do rings on the tools sit in the groove on the tool shank?

M Check all tools and push rings that have moved into the groove again.

 Do the parameters of the job in DentalCAM correspond to the parameters of the blank?

M Make sure that the following parameters of the job and of the blank correspond to each other and that they are suitable for the objects that you want to machine:

- Material type
- Dimensions of the blank
- Indications (types) of the individual objects

 Do you use the newest version of the manufacturing software that is released for the N4 Impression?

 New versions often contain important optimizations.

S  Update the applications by clicking on the depicted icon. Alternatively, contact customer service.

 Are the STL files of sufficient quality?

S1. Before you import STL files into the manufacturing software, check their quality. Especially consider the manufacturer information on wall thickness and border thickness.

S2. If necessary, adjust the settings of your scanner and scan application.

 Are the bars of the individual objects set correctly?

S If you set bars, consider the corresponding information in the documentation for DentalCAM.

 Is the measuring key contaminated? This prevents a correct tool measurement.

M Clean the measuring key with a brush.

 chapter 6.4, page 25

 Is the collet chuck of the spindle contaminated or is it loosely seated in the spindle?

M1. Clean the collet chuck with the provided spindle service set.

M2. When you insert the collet chuck into the spindle, make sure that it is properly seated.

 chapter 6.5, page 26

 Did you calibrate the machine shortly before the problems occurred or have there been any problems recently with the tool change?

M Contact customer service.

 Did you exchange the collet chuck at the recommended interval?

M Verify the recommended interval for exchanging the collet chuck in the maintenance table and exchange the collet chuck if necessary.

 chapter 6.2, page 24

 **DentalCNC interrupts the job execution and displays that the flow rate is too low**

 You can resolve this problem while the job execution is interrupted. DentalCNC continues the job as soon as the flow rate is sufficient again.

 Is there enough cooling liquid in the cooling liquid tank?

M If necessary, fill a sufficient amount of cooling liquid into the cooling liquid tank.

 chapter 4.5, page 14

 Is the pump defective?

M Replace the pump against an original manufacturer's spare part.

 chapter 6.10, page 30

 **DentalCNC interrupts the job execution and displays that the air pressure is too low.**

 You can resolve this problem while the job execution is interrupted. DentalCNC continues the job as soon as the air pressure is sufficient again.

 Is the service unit properly set?

M Set the air pressure at the service unit to a value between 4 bar and 8 bar (recommended: 4.5 bar).

 chapter 4.4.4, page 14

 Is the error caused by the external compressed air supply?

- M1. Close the external compressed air supply valve.
- M2. Check if all pneumatic hoses are properly seated in their connections and are undamaged.
- M3. Check if your compressor is switched on and set correctly.
- M4. Check if all effected valves of your compressed air system are open.

 Does the air pressure fluctuate significantly according to the service unit so that the job execution is interrupted and continued frequently?

- M1. Check if your compressor can *permanently* generate at least 4 bar air pressure with a volume flow of 35 l/min.

 **Not every compressor is designed for commercial use with dental machines.**

- M2. If necessary, replace your compressor with one that meets the above recommendations.

 chapter 4.4, page 13

 **I have exchanged the tool magazine inserts and now there are no holes for the tools anymore.**

 Tool magazine inserts are delivered without drilled holes for tools. You drill them with the machine.

- S Drill the tool positions into the inserts with DentalCNC.

 chapter 6.15, page 34

Index

A		
Abutment holders, using.....	20	
Air purity.....	13	
Axis calibration.....	32	
B		
Bellow.....	9	
Blank holder.....	9, 19	
Blanks, mount/unmount.....	19	
C		
Calibration.....	32	
Calibration specimen.....	33	
Cleaning.....	6, 24	
Collet chuck.....	9	
Clean.....	26	
Colours working chamber lighting.....	9	
Compressed air		
Check.....	16	
Consumption.....	13	
Install.....	13	
Requirements.....	13	
Set via service unit.....	14	
Condensate.....	31	
Connection panel.....	8	
Cooling liquid		
Exchange.....	29	
Fill.....	14	
Flow rate, check.....	16	
Cooling liquid container.....	9	
Clean.....	29	
Fill.....	14	
Remove from machine.....	28	
Cooling lubricant.....	10	
Mix.....	14	
D		
Discharging screw.....	31	
Disposal		
Cooling liquid.....	35	
Machine.....	35	
E		
Emergency release of the working chamber door.....	22	
F		
FAQ.....	36	
Filter cartridge.....	31	
Flow rate, check.....	16	
Fuse, exchange.....	31	
G		
General safety instructions.....	4	
H		
Housing, clean.....	30	
I		
Inspection.....	25	
Installation site.....	11	
Requirements.....	11	
Installing the machine.....	11	
Intended user.....	6	
ISO 8573-1.....	13	
J		
Job execution, start.....	21	
K		
Key for emergency release.....	22	
M		
Machine malfunction.....	22	
Main fuse, exchange.....	31	
Maintenance.....	6, 24	
Table.....	24	
Manufacturing software.....	10	
Measuring key.....	25	

Mixing ratio water emulsion 14

O

Operating regulations 6

Operation 17

Operation interruption 22

P

Power failure 22

Power outage 22

PreFace® abutment holder 20

PreFace® abutments 20

Pressure (compressed air)
Check 16, 37

Processing, start 21

Pump exchange 30

Pump filter 30

R

Room temperature 7, 11

Rotational axis 9, 10

S

Safety interlock 22

Scope of delivery 11

Service unit

Check 31

Install pneumatic hoses 14

Overview 13

Set air pressure 14

Sound emission 10

Spindle 7, 10

Storage 7

T

Technical data 10

Testing the machine 16

Test specimen 33

Tool breakage 22

Tool magazines 10, 18

Change inserts 34

Tool positions 18

Tools, chuck in/exchange 18

Transport 7

Troubleshooting 36

W

Warranty 24

Water emulsion

Exchange 29

Fill 14

Water separator 31

Working chamber

Clean 25

Lighting 9

Overview 9

Working chamber door, emergency release 22

EC Declaration of conformity

according to EC directive for machinery 2006/42/EC Annex II A

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fulfils all the relevant provisions of the following directives:

- 2006/42/EC Machinery directive
- 2014/30/EU EMC directive

References of the applied harmonised standards according to article 7 passage 2:

- EN 614-1:2006+A1:2009
- EN ISO 12100:2010
- EN 13128:2001+A2:2009
- EN ISO 13849-1:2015
- EN ISO 13849-2:2012
- EN 60204-1:2006
- EN 61326-1:2013
- EN 61326-2-1:2013
- EN 61000-3-2:2014
- EN 61000-3-3:2013

The manufacturer undertakes to transmit, in response to a reasoned request by the national authorities, relevant information on the machinery electronically. Person being established in the community, who is authorised to compile the technical file:

Dipl.-Ing. (FH) Frank Benzinger
Vorstandsvorsitzender / Chief Executive Officer (CEO)
vhf camufacture AG
Lettenstraße 10
D-72119 Ammerbuch

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